

You have
this idea to
ORGANIZE.

Here's how to do it.


**THINK IN
COLOR!**



Color Management Guide



Zirc Colors



The Ultimate Guide To Organization

You want to be organized and you
need practical ideas you can use.

We have a product line with
brilliant style designed to bring
[organization] and [color] to every room.

Ready?

Set?

Get Inspired!

How Do You Use a Color Code System?

[1] Select a color for each procedure you practice.

Example:

Crown & Bridge = Neon Orange
Endo = Teal
Composite = Neon Blue

[2] If you are a group practice select a color for each doctor and or operator. (For a solo practice skip this step.)

Example:

Dr. Smith = Neon Purple
Dr. Mark = Neon Pink

[3] Assess the number of procedures you do weekly (on average).

Example:

Crown & Bridge = 10 per week
Endo = 5 per week
Composite = 12 per week

[4] Determine quantity of trays you will need for each procedure. Take your configured numbers from step 3 and divide them by 3 (always round your numbers up).

Example:

Crown & Bridge = 4 Neon Orange Trays needed
Endo = 2 Teal Trays needed
Composite = 4 Neon Blue Trays needed

What Products Do I Need For an Effective Color Code System?

Set Up Trays

Our exclusive “lock-on” tray covers provide an excellent method of infection control and safety that allow you to do two things that other trays don't:

- [1] This system allows you to set-up your trays ahead of time for the day's procedures and “lock on” a cover to prevent cross-contamination.
- [2] Prevents any contaminated items from falling off the used tray during transportation from operator to central sterilization. The contaminated items stay enclosed until you are ready to process trays for sterilization.



Material Management

The next step in setting up your color code system is organizing your materials by procedure. This starts by removing your materials from your operator drawers and utilizing a more efficient alternative, the Zirc “Procedure Tub”. You will need one Tub per procedure (group practices may demand additional).

All Tubs are suggested to be stored in a common area such as “central sterilization”. Or they may be stored in the operator. Tubs are used as a way to transport items from your central sterilization to the operator and back again. Materials used during the procedure are delivered from the Tub to the doctor, keeping your procedure tray free for instrumentation.



At the end of the day, materials are wiped down and Tub is stored. Both Trays and Tubs can be stored in Zirc's Racks.

Call to request a Material Management Guide for more information on organization. [\[800-328-3899\]](tel:800-328-3899)



Instrument Management

You have determined the quantity of Trays & Tubs your practice needs, now you need to color code and organize your hand instruments, burs and endo files (we refer to these items as “Instrument Holders”. Zirc offers a variety of “Instrument Holders” that will manage these instruments during the sterilization process and at chair side.

How Many Do I Need?

Order the same number of “Instrument Holders” for each procedure as you did Trays.

Example:

Crown & Bridge: You designated Neon Orange as the procedure color. You have determined 4 Trays are needed; therefore, you will need 4 Neon Orange Instrument holders and Bur holders. For more information please visit www.zirc.com.

Hand Instrument Organizers



E-Z Jett Cassette



Compact Cassette



Compact Steri Container



Steri Container



Mini Matt

Bur Organizers



8 Hole Bur Blocks



14 Hole Bur Blocks



12 Hole Bur Guards



22 Hole Bur Guards

Endo Organizers



Endo Guards (hand)



Endo Guards (engine)



Endo Assist

E-Z ID Rings and Tape

The final step with your color code set-up is utilizing color code rings and tape. E-Z ID Rings are made to slide onto your instruments. E-Z ID Tape may be used to color code anything that a ring will not fit (i.e. larger instruments, bur holders, cabinet drawers, tubs, etc...). Both can be used in the following ways:




- [1] **Solo Practice:** Place color code “E-Z ID Rings” on all hand instruments in the corresponding procedure color

Example:

You would place Neon Orange rings on all of your Crown & Bridge Instruments.

- [2] **Group Practice:** Group practices would use two different E-Z ID Rings on each hand instrument. One color would signify the procedure (Neon Orange is for Crown & Bridge) and the other color would signify the doctor (Neon Purple is for Dr. Smith).
- [3] **Instrument Order:** The E-Z ID ring should be placed at a diagonal. This indicates to the assistant the order in which the doctor will be using the instruments.





Congratulations! Your procedures are now organized and color-coded.

