

Zirc Dental Products

Product Catalog

Helping dental teams be safe, happy and efficient.



WE HELP DENTAL TEAMS BE
HAPPY, SAFE, AND EFFICIENT



COLOR METHOD 3-16

Color Method Before & After	3
All-In-One Kits	4
Build Your Own Tub	5
Tub Accessories	6
Trays & Covers	7
Multi-Mod Racks	8
Cassettes & Instrument Holders	9-10
Steri-Containers	11
Bur & Endo Guards/Blocks	12-14
Rings & Tape	15
Consulting Services	16

ISOLATION 17-21

Mr. Thirsty® One-Step	17
Airway Armor®	18
Pink Petal™	19
Insti-Dam®	20
Lingua-Fix® & SE Cushion	21

MOUTH MIRRORS 22-26

Crystal HD® Mirrors	23-24
Mirror Magic® Anti-Fog Solution	25-26

OPERATORY/MISC 27-37

Denture & Retainer Boxes	27-28
Drawer Organizers	29
Steri-Soaker & CAD/CAM Block Organizer	30
Cabinet Trays & Lab Pans	31
E-Z Shelf & Accessories	32
E-Z Storage & Handi-Hopper	33
HVE + SE Accessories & Needle Cappers	34
Water Wise & Dispensers	35
Dispos-A-Screen Traps & Cuspidor Screens	36
Mighty Mixer Bowls, Spatula, and Organizers	37



You may feel organized,
But are you efficient?

We'll teach you critical organization protocols to help you clear the clutter!

colormethod

ZIRC.COM | COLORMETHOD.COM



it's
time to

RE ORGANIZE

your instruments and materials

...and we can help!

DETAILS ON PAGE 16

Small Changes, Big Results

Stay efficient with color-coded organization.

Zirc's Color Method creates a streamlined, stress-free day-to-day routine from setup to sterilization. Implementing this simple system adds a visual element to both instrument and material management yielding a workflow properly optimized for efficiency, leaving your team more confident and less stressed.



Composite Kits S

Store and organize composite materials with our convenient Composite Kits.



20Z489
Syringe Composite Kit
\$75

Syringe holder allows you to store up to 40 syringe tubes.

Includes:

- Syringe 20-Unit Stand
- 1 Long Tub Cup/Cover
- 1 Divided Slide Tray
- 1 Syringe Clip-On Holder
- Shade Label Sheet 2pk



20Z479
Capsule Composite Kit
\$63

Long Tub Cups provide a total of 20 compartments for capsule storage.

Includes:

- Tub Liner
- 4 Long Tub Cup/Cover
- 1 Divided Slide Tray
- 1 Capsule Clip-On Holder
- Shade Label Sheet 2pk



All-In-One Kits

Zirc's procedure tubs are designed to store and organize materials to keep your inventory under control.

20Z455

Complete Tub


\$83 Available in all colors

Includes:

- Procedure Tub
- Tub Liner
- Tub Dividers
- Single Tub Cup/Cover
- Double Tub Cup/Cover
- Divided Slide Tray
- Safe-Lok® Cover for Tub

Dimensions: 12-5/8" x 11-1/4" x 4"

S




ZIRC TIP

Improve your clinical workflow!



The extensive variety of materials needed adds complexity to the restocking flow, making it challenging to streamline storage. By procedure, and in the same color as your trays, tubs now act as portable drawers for 'grab-and-go' efficiency.

- Syringe Composite Tub
- Capsule Composite Tub
- Endodontics Tub
- Crown Prep Tub
- Hygiene Tub
- Oral Surgery Tub






Build Your Own Tub





1 20Z463
Procedure Tub
\$22.50 Available in all colors
 Allows you to "build" your own setup.
 Dimensions: 12-1/4" x 10-7/8" x 2-3/4"
 

2 20Z462
Tub Liner
\$6.50
 Foam liner that sits in bottom of Tub.
 Tub not included. Replace as needed.
 Dimensions: 9-5/8" x 8-7/8" x 1/8"

3 20Z461
Tub Dividers
\$15 Gray only
 Creates up to 24 compartments. May be cut to create larger compartments.
 Dimensions: 9-5/8" x 8-5/8" x 1-5/16"


4 20Z459
Safe-Lok® Tub Cover
\$22.50
 Locking cover for securing materials. Keeps materials protected during storage and transportation.
 Dimensions: 12-5-8" x 11-1/4" x 1-7/16"
 

5 20Z465
Tub Insert
\$17.50 Gray only
 10 compartments for easy organization. Sits on rim of Tub, allowing for storage underneath.
 Dimensions: 10-3/4" x 9-11/16" x 1"


6 20Z474
Tub Box
\$16.50 Gray only
 Designed to hold/organize supplies. Non-slip bottom for secure placement. Flip-top cover.
 Dimensions: 3-3/4" x 3-5/8" x 1-7/8"


Tub Accessories

Tub Cups with Covers S

Fill Tub Cups with brushes, clamps, cotton pellets, prophylaxis angle tips, prophylaxis paste, pumice, syringe tips, implant tools, etc. Designed to be used inside the Procedure Tub.



20Z473
Long Tub Cup w/ Cover
\$15.50 Gray only

Internal dividers are removable.
Dimensions: 8-3/8" x 2-5/16" x 1-1/8"




20Z472
Double Tub Cup w/ Cover
\$9.50 Gray only

Dimensions: 3-3/8" x 2-1/4" x 1-1/8"



20Z471
Single Tub Cup w/ Cover
\$9.50 Gray only

Dimensions: 2-1/4" x 1-5/8" x 1-1/8"



20Z477
Shade Label Sheets 2pk
\$6

Preprinted labels for marking composite materials within the Tub organizers.

Dimensions: 7" x 4-3/8" x 1/64"

Slide Trays S ⌵

Sits on rim of Tub and allows for storage underneath.



20Z466
Divided Slide Tray
\$14.50 Clear


Dimensions: 9-3/4" x 5-3/8" x 1-1/16"



20Z467
Flat Slide Tray
\$14.50 Clear

Dimensions: 9-3/4" x 5-3/8" x 1-1/16"


Syringe Holders S ⌵



20Z487
Syringe 20-Unit Stand
\$26

Stands up for easy access to syringe-able materials. Lays flat to fit inside Tub or drawer. Holds up to 20 syringes. Has space for 2 Syringe 10-Unit Add-Ons (Up to 40 syringes).

Dimensions: 8-3/4" x 4-3/4" x 1-9/16"



20Z486
Syringe Clip-On Holder
\$11

Allows quick access to composite during procedure. Conveniently clips on any side of the Tub or B-Lok Tray. Holds 15 syringes and 3 capsules.

Dimensions: 8" x 2-7/16" x 3/4"



20Z488
Syringe 10-Unit Add-On
\$11

An add-on for the Syringe 20-Unit Stand. Add for storage of up to 40 syringe tubes. Gives additional space for syringe storage. Holds 10 syringes.

Dimensions: 8-9/16" x 2-5/16" x 3/4"



20Z476
Capsule Clip-On Holder
\$11

Allows quick access to composite during procedure. Conveniently clips on any side of the Tub or B-Lok Tray. Holds 8 capsules, and 4 syringes.

Dimensions: 8" x 2" x 9/16"

Trays & Covers

Our exclusive B-Lok trays and covers provide an excellent method of infection control and safety that allows for two benefits other trays do not: preventing any contaminated items from falling off the used tray during transportation from operator to central sterilization, and keeping contaminated items enclosed until you are ready for the cleaning and sterilization process.



20Z401
B-Lok Flat Tray
\$15 Available in all colors

Wide edges for secure grip. Non-slip feet for stability. High-gloss surface helps prevent staining.


Dimensions: 13-3/8" x 9-5/8" x 7/8"



20Z451
B-Lok Divided Tray
\$17 Available in all colors

Holds 12 hand instruments, including 9 FG, RA, or HP bur holders. Features Insti-Grip, non-slip feet for stability, and self-lubricating, stain-resistant surface. Also includes holder for Syringe Tip or Saliva Ejector.

Dimensions: 13-3/8" x 9-5/8" x 7/8"



20Z454I
Insti-Grip
\$3

To be used with the B-Lok divided Tray only. Holds instruments quietly and safely in place.

Dimensions: 6-5/8" x 1/4" x 3/8"



20Z101
Mini Tray
\$13.50 Available in all colors

Wide edges for secure grip. Non-slip feet for stability. High-gloss surface helps prevent staining. Mini Trays are great for hygiene/exam setup.

Dimensions: 9-3/8" x 6-3/8" x 7/8"



20Z445
Safe-Lok® Tray Cover
B-Size
\$19.50

Fits B-Lok Flat Tray and Divided Trays only. Slide locking mechanism securely locks cover to tray.

Dimensions: 13-7/8" x 9-7/8" x 3/4"



Tall B-Size Standard B-Size

20Z446
Safe-Lok® Tray Cover
Tall B-Size
\$20.50

The Tall size accommodates our Performance Tool Organizer and almost any stainless steel cassette. Fits B-Lok Flat and Divided Trays only. Side locking mechanism securely locks cover to tray.

Dimensions: 13-7/8" x 9-7/8" x 1-3/8"



20Z441
Tray Cover B-Size
(Non-Locking)
\$19

Dimensions: 14" x 10-1/8" x 5/16"



20Z103
Safe-Lok® Mini Tray Cover
F-Size
\$17.50

Side locking mechanism securely locks cover to tray.

Dimensions: 9-7/8" x 6-5/8" x 1"




20Z102
Mini Tray Cover F-Size
(Non-Locking)
\$17

Dimensions: 9-5/8" x 6-5/8" x 7/8"



Multi-Mod Racks

Tub and tray racks are effective additions to a dental practice. These racks have non-slip feet so they are stable on any counter or cabinet surface. With a chrome-plated finish, these racks can fit any dental practice design.



21Z410
Multi-Mod Rack 8-Place
\$155

Holds 8 Trays or 4 Tubs with covers. Sturdy metal construction; chrome-plated. Non-slip feet. Adjustable to hold trays and tubs.



21Z105
Multi-Mod Rack 6-Place
\$148

Holds 6 trays or 3 tubs with covers. Adjustable to hold trays and tubs. Sturdy metal construction; chrome-plated. Non-slip feet.



21Z101
Mini Tray Rack
\$120

Holds 6 mini trays with covers. Chrome-plated, sturdy metal construction with non-slip feet. Fits size F trays.



ZIRC TIP

Are you in compliance?

The Safe-Lok® Tray Cover helps in compliance with OSHA and CDC,

OSHA

Immediately or as soon as possible after use, contaminated reusable sharp shall be placed in a puncture-resistant, leak-proof and labeled or color-coded containers until properly reprocessed. To avoid spillage of contents, it is suggested that they be covered and secured prior to moving.

CDC

Contaminated instruments must be "contained" prior to transporting them. Containment can be defined as instrument cassettes or cages in which dental care items are cleaned, sterilized, and stored until point of use. Sharp instruments should not be carried openly to the instrument-processing area. Percutaneous injuries may occur to other dental healthcare personnel (DHCP) or patients in hallways leading to the instrument-processing area.



Why Use Zirc's Cassettes?

A low-cost investment with a quick return! Individual instrument handling takes time away from treating patients and increases overall reprocessing time. Zirc's cassettes improves instrument efficiency, safety, and compliance without breaking the bank.

Product Benefits

- Zirc's cassettes will save you time and help prevent injury.
- Keep instrument setups together; no sorting.
- Protects your instruments from damage and prolongs their life.
- Prevents the team from a potential exposure.
- Better patient perception, making patients less anxious.

Product Features


- Designed for superior water flow and cleaning – open top, round holes, and radius corners.
- Keeps instrument setups together through sterilization and procedure processes.
- Non-slip feet for stability and safety.
- Autoclavable
- Available in 16 colors for color coding by procedure.

Sterilization Pouch Size

Cassette	Pouch Size
E-Z Jett® 5-Place	6" x 12"
E-Z Jett® 8-Place	7 ½" x 13"
E-Z Jett® 8-Place Slim	6" x 12"
E-Z Jett® 10-Place	10" x 13" or 14"
Performance Tool Organizer	7 ½" x 13"




Cassettes & Instrument Holders



50Z930
E-Z Jett® Cassette
10-Place
\$49 Available in all colors

Holds 10 instruments. Accessory compartment in the middle for extra storage. Covers stay attached to the cassette and fold back easily.


Dimensions: 7-15/16" x 7" x 1-1/8"



50Z927
E-Z Jett® Cassette
8-Place
\$44 Available in all colors

Holds 8 instruments. Covers stay attached to the cassette and fold back easily. A/W Syringe Tip Holder included.

Dimensions: 7-15/16" x 4-7/16" x 1-1/8"




50Z911
E-Z Jett® Cassette
Slim 8-Place
\$39 Available in all colors

Holds 8 instruments. Covers stay attached to the cassette and fold back easily. Compact size. Fits the Statim. A/W Syringe Tip Holder included.

Dimensions: 6-7/8" x 3-9/16" x 1-3/16"


* Size comparison with the 8-Place E-Z Jett® Cassette



50Z925
E-Z Jett® Cassette
5-Place
\$34 Available in all colors

Holds 5 instruments. Great for basic exam setup! A/W Syringe Tip Holder included. Covers stay attached to the cassette and fold back easily. Fits Statim Sterilizer.


Dimensions: 7-5/16" x 3-1/8" x 1-3/16"



20Z464
Instrument Mat
\$16.50 Available in all colors

Holds 10 hand instruments. Reverse side holds 4 surgical instruments. Non-slip, soft material. Organizes instruments in order of procedure use. Can be cut to fit Mini Tray.

Dimensions: 6-3/4" x 2" x 3/8"



20Z480
Mini Mat
\$15 Available in all colors

Holds 8 instruments. Organizes instruments at chairside and for sterilization. Non-slip feet for stability.

Dimensions: 3-7/8" x 3-1/2" x 1/2"


50Z917
Instrument Strap
\$6.50/Pair Gray only

Locking strap secures instruments.

Dimensions: 3-5/8" x 1/2" x 5/16"

50Z926
Syringe Tip Holder 5pk
\$5.25

Dimensions: 1" x 3/8" x 3/8"



Performance Tool Organizer



A deeper, roomier cassette made for large and specialty instruments. Designed for superior water flow for superb cleaning. Compact styling to fit in most sterilizers. Covers stay attached to the cassette and fold back easily. Comes with two strips of 10 customizable removable feet to design your own layout.

50Z940

Performance Tool Organizer

\$63 Available in all colors

Dimensions: 8-3/16" x 4-1/4" x 1-7/32"

50Z944

Performance Tool Organizer Replacement Feet

\$8 Gray only | Two strips of 10.



50Z900

Steri-Container

\$13 Available in all colors

Keeps instrument setups together through sterilization processes. Holds 10-14 instruments.

Dimensions: 8-1/8" x 1-7/8" x 1-7/8"



50Z905

Compact Steri-Container


\$13 Available in all colors

Easy open, push down latch. Keeps instrument setups together through sterilization processes. Statim 2000 holds 7 Compact Steri-Containers. Holds 8-12 instruments.

Dimensions: 7-1/8" x 1-1/2" x 1-1/2"



ZIRC TIP



Closed

Open

Hold your large & specialty tools!

Performance Tool Organizer fits on Zirc Trays with the Safe-Lok® Tall Cover!

- Ideal for:
- Surgical elevators and forceps
 - Orthodontic pliers
 - Hemostats
 - Scissors
 - Any instrument that does not fit in your regular cassette



Bur Blocks

Our Bur Blocks will organize and keep your burs safe during the sterilization process.

Magnetic Bur Blocks

- Fits FG, RA, and HP Burs
- Magnetic base holds burs securely
- Non-slip feet for stability
- Drain holes help prevent rusting of burs



50Z400
Magnetic Bur Block 8-Hole
\$16.50 Available in all colors

Covers sold separately.

Dimensions: 1-1/2" x 1-1/2" x 1/2"



50Z403
Magnetic Bur Block 14-Hole
\$17 Available in all colors

Covers sold separately.

Dimensions: 2-15/16" x 1-1/2" x 1/2"



Recessed area to hold used burs.



Bur Block Covers

- Covers available to maintain cleanliness of burs during storage
- 1" covers for FG and RA burs; 2" cover (14-hole only) for lab and surgical burs
- Covers are NOT autoclavable
- Surface disinfect only



50Z401
Magnetic Bur Block 8-Hole Cover 1" High Standard
\$8.50

Dimensions: 1-5/8" x 1-5/8" x 1-1/16"



50Z404
Magnetic Bur Block 14-Hole Cover 1" High Standard
\$8.50

Dimensions: 3-1/8" x 1-5/8" x 1-1/16"



50Z405
Magnetic Bur Block 14-Hole Cover 2" High HP Lab Burs
\$8.50



Dimensions: 3-1/8" x 1-5/8" x 2-1/4"

Steri-Bur Guards & Adapters

This gadget has it all! Not only does it protect your burs during cleaning and sterilization, but it can also be positioned at three different angles to give you easy access at chairside.

Steri-Bur Guard

- Fits FG, RA, and HP Burs.
- Adjustable height for burs of different lengths (12-Hole Tall Bur Guard not adjustable).
- Holds burs in position for maximum cleaning effectiveness.
- Drain holes help prevent rusting of burs.



50Z406
Steri-Bur Guard 12-Hole

\$17 Available in all colors

Accommodates burs from 12mm-30mm.

Dimensions: 2-7/8" x 3/8" x 1-3/8"



50Z410
Steri-Bur Guard 22-Hole

\$18 Available in all colors

Accommodates burs from 12mm-30mm.

Dimensions: 4-7/8" x 7/16" x 1-1/2"

50Z412
Universal Bur Adapter 5pk

\$5.25
Can be used with 12 or 22-Hole Bur Guard. Enables long and short burs to be used in the same bur guard without short burs falling out.

A short bur adapter is included to stop any shorter burs from falling out.

50Z408
Steri-Bur Guard Tall 12-Hole

\$17.50 Gray only

Accommodates burs from 30.9mm-37.1mm.

Dimensions: 2-7/8" x 3/8" x 1-5/8"



Steri-Endo Guards & Holders

A unique delivery system to organize, store, and sterilize endodontic files. The Endo Guard swivels 360° with stops at every 45° for operator's choice. Large slotting provides complete visibility and draining during cleaning.

Steri-Endo Guard

- Holds most lengths for staging, cleaning, and sterilizing files and reamers.
- Millimeter gauge measures and positions rubber stopper.
- Counter on locking bar tracks the number of sterilizations.
- Right- or left-hand use; optional angling.
- Available in Hand or Engine styles.
- Single-use foam insert cleans/holds files (packaged with 12).



50Z450
Steri-Endo Guard Hand
\$35.50 Available in all colors
 Organizes up to 16 hand files/reamers.
 Dimensions: 5-3/8" x 3/8" x 2-3/16"

50Z451
Foam Inserts 48pk
\$35.50 Single use

50Z451V
Foam Inserts 144pk
\$70 Single use



50Z452
Steri-Endo Guard Engine
\$35.50 Available in all colors
 Organizes up to 11 engine and 5 hand files/reamers.
 Dimensions: 5-3/8" x 3/8" x 2-3/16"



50Z460
Endo Assist
\$31.50 Available in all colors

The Endo Assist features a built-in millimeter gauge for precise measuring and positioning of rubber stoppers. Designed for both right and left-handed use, it includes small and large rings to fit most finger sizes. Its flat bottom allows it to stand independently, and it comes with a single-use foam insert for file cleaning and storage (packaged w/ 12).

Dimensions: 1-1/2" x 1-3/4" x 1-5/8"



50Z462
Endo Assist Stand
\$18 Black only

Holds endo files and reamers in upright position. Foam provides fast and easy cleaning of the files during use.

Dimensions: 1-7/8" x 1-7/8" x 1-1/4"



50Z466
Endo Assist Stand Foam 48pk
\$34.50 Single use

Dimensions: 10-1/2" x 7-3/4" x 7/8"

50Z466V
Endo Assist Stand Foam 144pk
\$70 Single use

Dimensions: 10-1/2" x 7-3/4" x 7/8"



Rings & Tape

How to effectively identify your instrument setup

1 Color for Procedure

Vibrant Orange is for identifying the procedure

1 Color for Practitioner (optional)

Vibrant Blue shows ownership to the practitioner

ZIRC TIP

Place your instrument to determine which ring size you need.

SM
3.17mm-6.35mm


LG
6.35mm-9.52mm

XL
9.52mm-12.7mm

Did you know?
The diagonal setup of the orange rings shows the order of how the instruments should be used. An optional second color may be used to identify the practitioner.

E-Z ID Rings

- Prepackaged in an easy-to-use dispenser
- Available in Small, Large, and Extra Large sizes to fit all instrument handles



Single Color 25pk

Available in all colors

- SM** \$8 | 70Z100
- LG** \$8.25 | 70Z200
- XL** \$8.50 | 70Z000

Classic System 200pk

Assorted colors [A, B, D, G, I]

- SM** \$52 | 70Z105
- LG** \$53 | 70Z205
- XL** \$54 | 70Z005

Vibrant System 200pk

Assorted colors [N, O, P, R, Q, S]

- SM** \$52 | 70Z110
- LG** \$53 | 70Z210
- XL** \$54 | 70Z010

Jewel System 200pk

Assorted colors [E, J, T, I, M, U]

- SM** \$52 | 70Z115
- LG** \$53 | 70Z215
- XL** \$54 | 70Z015

E-Z ID Tape

- Tape dispenser has built-in tape cutter
- Dispenser roller prevents tape from twisting



E-Z ID Tape System

\$49 Includes 8 - 3' roll

Classic | 70Z305

Assorted colors [A, B, D, G, I]

Jewel | 70Z315

Assorted colors [E, J, T, I, M, U]

Vibrant | 70Z310

Assorted colors [N, O, P, Q, R, S]

70Z300

E-Z ID Tape Roll

\$12 Available in all colors

Includes 1 - 10' roll.



Consulting Services

Color Method will rescue your team from clinical clutter and disorganization. It's time to introduce a color-coded organization system to your practice. Schedule a consultation or reset and we'll do it for you!

BRONZE General Guidance

This package is perfect for any dental professional to get a quick win and learn more about lean clinical performance. We'll also explain our process, recommend a next step, and answer your top questions.

- 25-minute consultation with a Clinical Efficiency Coach
- Learn the secret to clinical organization
- Get your questions answered

Free

SILVER Clinical Efficiency Consultation

This package is for anyone ready to improve overall clinical performance. We'll send you a demo box full of our products to experience for yourself. We'll take the time to get to know you and your space and then present you with actionable steps to take on your own.

- Professional clinical organizer for 4-6 hours
- Custom materials playbook
- Demo kit to explore products
- 15% off first order
- Custom plan for lean clinical performance
- 1-hr. team training
- Team members earn 1 CE credit

\$647

GOLD Virtual Color Method

You've decided you absolutely need Color Method, the whole team is on board and you're ready to achieve near-perfect clinical performance. In this package, we will design the Color Method system for you, help you build your tubs and trays and be a part of the Implementation Day to guide you to success!

- Professional clinical organizer(s) for 10-12 hours
- Custom plan for lean clinical performance
- Full custom playbook
- 1-hr team training
- Demo kit to explore products
- Team members earn 3 CE credit

\$1900

PLATINUM On-Site Color Method

Time is low. Your team is stressed. There is no margin for error and the practice needs to stay open as long as possible. It's called Platinum for a reason: our clinical consultants specialize in meticulously tailoring every aspect of your process, and they'll be there with you in person on Implementation Day to provide hands-on coaching.

- Professional clinical organizer(s) for 18-24 hours
- Custom plan for lean clinical performance
- Full custom playbook
- 1-hr team training
- Demo kit to explore products
- Review OSHA compliance
- Team members earn 3 CE credit
- 15% off orders for 1 year

\$4900

QUARTZ Full Practice Evaluation

This package is ideal for any denovo or renovation project. We'll review your floor plans, technology, and storage areas to minimize mistakes and regrets for your build. Our recommendations will save hours of time with your architects and contractors, and pave the way for lean clinical workflows.

- Complete clinical walkthrough and evaluation
- Deep assessment of all cabinetry and storage areas
- Instruments & materials analysis
- Review OSHA compliance
- Team members earn 3 CE credit
- 15% off orders for 1 year

Contact For Pricing

DIAMOND Platinum & Quartz

For a nice discount, combine Platinum and Quartz and know with absolute confidence that your clinical processes will become effective, leading to a more efficient and harmonious clinic.

- Hit reset on your entire clinical operations.
- Combine our full service, on-site clinical consultation services and transform your practice
- Team members earn 3 CE credit
- 15% off orders for 1 year

Contact For Pricing

Need help deciding? Let us help you!


Call (763) 682-6636 or visit www.colormethod.com



mr. thirsty®

one - step

Hands-free retraction, isolation, and evacuation provided by one simple device! Attach Mr. Thirsty® to your existing HVE, perform the procedure, and dispose when done! Rather than time spent on sterilizing or tracking number of uses, it can be spent on additional procedures, daily tasks, or catching a break! According to a time and motion study, Mr. Thirsty® can reduce chair procedure time by up to 29%!



50Z992
Mr. Thirsty® Comfort Kit
\$65

The Comfort Kit's lightweight tubing allows your heavy HVE Valve/Tubing to stay mounted to the equipment.

Includes:
 Lightweight lever HVE Valve,
 6' lightweight tubing.

Mr. Thirsty® One-Step NOT included.



Trimmable

All three sizes can be trimmed to fit, eliminating the need for overstock



Comfortable

Offers tongue and cheek retraction with a built-in bite block for comfort.



Disposable

No time wasted in maintenance or sterilization



Adaptable

Connects to existing HVE, allowing it to be used in any operator



Affordable

No equipment investment

Mr. Thirsty® One-Step

- Keeps work environment dry
- Fully assembled
- Single use
- Latex-free

50pk
\$125

- 50Z989SM | Small
- 50Z989MD | Medium
- 50Z989LG | Large

100pk
\$225

- 50Z989SM-B | Small
- 50Z989MD-B | Medium
- 50Z989LG-B | Large

50Z988
Variety Pack
\$84 Includes 10 of each size

- (10) Small
- (10) Medium
- (10) Large



Smaller size is perfect for smaller mouths!





AIRWAY ARMOR®

With all the small items used in dentistry, (burs, files, implants and veneers) teams always run the risk of something being aspirated or ingested. Airway Armor® was created to eliminate the unnecessary risk for patients and protection for the practice.

Ideal for Crown & Bridge, Endodontics, Extractions, Implants, Inlays/Onlays, Orthodontics, Veneers

- Aspiration and ingestion prevention
- Patient protection and comfort
- Easy placement and removal
- 3 convenient sizes – small, medium, and large
- Single use
- Latex-free

Airway Armor® 25pk

\$25

50Z960SM | Small

50Z960MD | Medium

50Z960LG | Large



50Z960TRI


15pk Trial Kit

\$20 Includes 5 of each size

(5) Small

(5) Medium


(5) Large



Easy placement



Protection



Pink Petal™



comfortable,
hands-free suction

Created by hygienist Karrie Velky, the Pink Petal™ is a great product for all dental practices looking to fulfill the need of hands-free, low volume evacuation. Easily slide the Pink Petal™ on any saliva ejector, creating a device that will provide continuous suction.

A wonderful adjunct tool to free up a clinician's hand

- Holds low-volume suction tip in place
- Protects and retracts the cheek
- Patient can close around suction with ease
- Easily slides onto saliva ejector
- Soft, comfortable material
- Single use
- Latex-free



50Z975
Pink Petal™ 50pk
\$27

INSTI-DAM®

Features for BOTH Relaxed Fit and Latex-Free Insti-Dam®

- Built-in flexible medium gauge material that comes pre-assembled with frame and pre-punched hole
- Pre-punched hole makes positioning easy and allows eyes and nose to be open and comfortable
- Additional holes can be punched
- Compact design fits outside patients' lips, easy-to-use, non-threatening, and comfortable for patients
- Radiographs may be taken by folding it to the side. No need to remove for mid procedure x-ray
- Single use
- Latex-free
- 50% faster placement than traditional dam.



50Z457
**Latex-Free Insti-Dam®
Relaxed Fit 20pk**
\$76

Ideal for posterior isolation. The extra "relaxed" material in the Relaxed Fit Insti-Dam® makes for quick posterior placement – no stretching, tearing, or pulling.

Dimensions: 4-7/8" x 4-7/8" x 6-1/8" diameter

Patent #7044738



50Z459
Latex-Free Insti-Dam® 20pk
\$66

Ideal for anterior isolation. Made with polyisoprene.


Dimensions: 4-1/4" diameter



50Z471
Insti-Dam® Dispenser
\$21 Only available in white

Can be mounted to wall or cabinet or can sit on a counter. Non-slip rubber bottom for stability. Holds 35 Insti-Dams®. Made for Latex-Free Insti-Dam® ONLY.

Dimensions: 4-7/8" x 4-7/8" x 6-1/8"



Lingua-Fix®

- Lingua-Fix® provides a dry work area
- Holds tongue in a steady position to protect from injuries
- Simple operation – connects to low-volume suction valve
- Comfortable, no sharp edges or corners
- Large suction surface avoids blockage and is more comfortable
- Reduces sedimentation of drilling byproducts in pipes and tubing
- Single use
- Comes pre-assembled




50Z953S
Lingua-Fix® 10pk
 \$14.50

50Z953
Lingua-Fix® 50pk
 \$51

SE Cushion

- Used with a Saliva Ejector
- Prevents tissue grabbing and irritation
- Increases patient comfort
- Quiets the noise of suction
- Single use

50Z970
SE Cushion 100pk
 \$24



Helpful Tip

Ideal for use with rubber dams

Comfortable and efficient suction or evacuation



CRYSTALHD®

What a view!

The brighter the workspace, the better the illumination and visibility. With Zirc's Crystal HD® lens, you are visually working in a brighter environment.

If you haven't made the change to HD – what are you waiting for?

As a microscope user, I can understand the importance of clear vision. Zirc mirrors significantly improved the brightness of my image, giving my eyes highest quality vision. It's true – brightness makes all the difference.


Dr. Dimitris Tsanaktsidis
Microscopic & Aesthetic Dentistry



Brighter than rhodium mirrors!
Image provided by Dr. Lori Trost


Why Use Zirc's Resin Mirrors?

Crystal HD® reflects color 99% accurately while Rhodium only reflects color 75% accurately.



	Resin	Stainless
Autoclavable	X	X
Lens flush with mirror frame	X	
No galvanic shock on patient	X	
Rust-free	X	

CUTTING EDGE IMAGERY



ZIRC TIP

What's Your Mirror Size?



diagram is to scale

- **#0 Mirror**
14mm
- **#3 Mirror**
20mm
- **#4 Mirror**
22mm
- **#5 Mirror**
24mm

Crystal HD® Mouth Mirrors



Crystal HD® Mouth Mirror Soft Grip 12pk
\$84

Assorted Vibrant
2 of each color
[N, O, P, Q, R, S]

#4 | 50Z363-NEON
#5 | 50Z364-NEON

Assorted Jewel
2 of each color
[E, J, T, I, M, U]

#4 | 50Z363-JWL
#5 | 50Z364-JWL



Crystal HD® Mouth Mirror Soft Grip Black/Gray

\$31 4pk \$84 12pk
#4 | 50Z361 #4 | 50Z363
#5 | 50Z362 #5 | 50Z364



Crystal HD® Mouth Mirror Thin Grip 12pk
\$84

Assorted Vibrant
2 of each color
[N, O, P, Q, R, S]

#4 | 50Z353-NEON
#5 | 50Z354-NEON

Assorted Jewel
2 of each color
[E, J, T, I, M, U]

#4 | 50Z353-JWL
#5 | 50Z354-JWL

MEGA Mirrors

*NOT a medical device

50Z309

MEGA Mirror

\$21 Color code: N, P, Q, R, S

A fun toy/prop for your patients.
2-1/2" lens and 15-1/2" in length.

S



Crystal HD® Mirror Heads  



Crystal HD® Resin Mouth Mirror Screw-on Heads - Cone Socket 12pk
\$84

Single Color
Available in all colors and K (black)
#4 | 50Z371
#5 | 50Z372

Assorted Vibrant
2 of each color:
[N, O, P, Q, R, S]
#4 | 50Z371-NEON
#5 | 50Z372-NEON

Assorted Jewel
2 of each color:
[E, J, T, I, M, U]
#4 | 50Z371-JWL
#5 | 50Z372-JWL

Assorted Classic
2 of each color:
[A, B, D, G, I, K]
#4 | 50Z371-PSTL
#5 | 50Z372-PSTL



Crystal HD® Resin Mouth Mirror Double-Sided Heads - Cone Socket 4pk
\$45

Assorted Vibrant
1 of each color
[N, O, R, S]
#4 | 50Z374-NEON

Gray Mirror Head
#4 | 50Z374



Crystal HD® Stainless Mouth Mirror Screw-On Head - Cone Socket
\$84

Stainless Steel 12pk
#3 | 50R379
#4 | 50R380
#5 | 50R381



Stainless Steel Head
\$14 Each
#0 | 50R378




Crystal HD® Mirror Handles  



50Z340
Soft Grip Handle - Cone Socket
\$18 Each Available in all colors



50Z342
Ergo Grip Handle - Cone socket
\$18 Each Available in K (black) & all colors except A



Mirror Magic® Anti-Fog System

The Mirror Magic® Anti-Fog System prevents fogging on mouth mirrors and offers the practitioner a quick method for removing debris during procedures. The disposable swipe pads stick to the back of a gloved hand for easy access. The Anti-Fog System requires only a small amount of solution that provides 45 minutes of fog-free working time and is economic in size, making it easy to store.



The Mirror Magic® Anti-Fog System decreases fogging while also decreasing the build-up of debris on my mirror. This affords me outstanding visualization, which is crucial to performing quality dentistry. mirror magic® has truly made a difference in my practice and is one of those features I didn't realize I needed; and now I can't imagine practicing without it.

Andrew M. Goldsmith
DDS, DICOI, FIALD

50R200
Mirror Magic® Trial Kit
\$26

Comes with 2-6ml bottles of Anti-Fog Solution and 24 Swipe Pads.



Without mirror magic®



With mirror magic®



The Perfect Pair

The superior visual acuity and brighter Crystal HD® lens combined with the quick debris removal and Anti-Fog Solution of mirror magic® have proven to be the perfect pair. Increase chairside efficiency and speed up your process with this dynamic duo.



50R201
Mirror Magic® Organizer
\$16.50

Holds: 48 Swipe Pads and 1-6ml bottle of Anti-Fog Solution. Store on operatory counter or in drawer. Non-slip bottom for secure placement.

Dimensions: 3-3/4" x 3-5/8" x 1-7/8"

Note: box only, does not include Mirror Magic®




50R210
**Mirror Magic®
Swipe pads 96pk**
\$26 Single use



50R205
**Mirror Magic®
Anti-Fog Solution 10pk**
\$76

Each bottle contains 6ml each and lasts approx. 50 patients.



Denture & Retainer Boxes

Retainer & Denture Boxes S

One piece construction. Solid locking mechanism.

Retainer Box 12pk

1" Deep | 25R500

\$14

Dimensions: 3-1/8" x 2-3/4" x 1"

1 1/2" Deep | 25R550

\$15

Dimensions: 3-1/8" x 2-3/4" x 1-1/2"

Single Colors



Assorted Pack Contains 3 of each color below



30R800

Denture Box 12pk

\$16

Dimensions: 4-3/8" x 3-3/8" x 2"

Single -or- Assorted Pack

Contains 3 of each color below



Retainer boxes can be hooked to backpacks for better portability!



Personalized Imprinting

Our new and improved imprinting service is quick and easy. Upload your own logo or select a full color icon from our pre-made Fun or Professional sets. Add up to four lines of text for your business name, phone number, website, or your own unique message. It's as easy as that!

Imprinted Retainer Box

24pk

1" Deep | 25P500

\$38

Dimensions: 3-1/8" x 2-3/4" x 1"

1 1/2" Deep | 25P550

\$40

Dimensions: 3-1/8" x 2-3/4" x 1-1/2"

1000pk

1" Deep | 25PB500

\$1,351

Dimensions: 3-1/8" x 2-3/4" x 1"

1 1/2" Deep | 25PB550

\$1,422

Dimensions: 3-1/8" x 2-3/4" x 1-1/2"

Single Colors



Assorted Pack Contains 6 of each for 24pk & 250 of each for 1000pk



- 1 Effortlessly upload your logo or choose from pre-made icons to imprint with our high quality, 4-color printing process.
- 2 Personalize with up to 4 lines of text
- 3 Manufactured and imprinted in the USA!



Imprinted Denture Box

24PK

30P800

\$42

Dimensions: 4-3/8" x 3-3/8" x 2"

1000PK

30PB800


\$1,460

Dimensions: 4-3/8" x 3-3/8" x 2"

Single Colors



Assorted Pack Contains 6 of each for 24pk & 250 of each for 1000pk



Drawer Organizers S

Zirc is excited to announce a new product – a modular Drawer Organizer. All Zirc products are designed with time in mind for an organized and efficient dental practice, this new Drawer Organizer is no different.

20Z210

Drawer Organizer


\$38 Color Code: A, G, I, N, P, R, S

Each Drawer Organizer will come with 1 large divider, 2 medium dividers and 4 small dividers. Additional dividers can be purchased.

Dimensions: 11.5" x 6.6" x 2"

Large Divider (Included)

Dimensions: 0.95" x 1.97" x 11.35"



Drawer Organizer dividers

\$10 Color Code: A, G, I, N, P, R, S

Dividers may be ordered separate. Size (S) comes in 2 pieces and size (M) in 1 piece

20Z207 | Small

Dimensions: 0.63" x 1.96" x 3.19"

20Z208 | Medium

Dimensions: 0.63" x 1.97" x 6.46"



The Steri-Soaker

An efficient, high-volume germicide tray designed for maximum use of chemicals to clean/disinfect instruments. It is designed to hold instruments for easy cleaning and is large enough to accommodate Zirc Steri-Containers, Cassettes, Bur Blocks, and Bur Guards.


50Z805
STERI-SOAKER

\$105 Color code: B, I

Calendar on lid conveniently reminds staff when to change the solution; use color code rings to mark dates. Transfer basket's ribbed and perforated design allows quick drain-off and easy pick-up of instruments. Holds 3.78L of liquid; fill lines shown in quarts and liters. Hinged fitted cover.

Dimensions: 13-7/8" x 9-5/8" x 5-1/8"

Safe to use chemicals: cleaners (enzymes, detergents), disinfectants (phenols)/cold sterilants for glutaraldehyde and chlorines.




50Z417
CAD/CAM Block Organizer
\$55 Color Code: E, J, N, Q, S, T, U

Provides a professional platform to present and discuss same-day restoration. Fits all CAD/CAM milling blocks including PrograMill®, CEREC®, PlanMill™ Blocks. Store in drawer or cabinet.

Dimensions: 5-9/16" W x 1-7/8" x 7-3/4" D



Pick your color!



Cabinet Trays & Lab Pans

Cabinet Trays



20Z201
**No 16A
Hand Instruments
(Regular Well)**

\$9.50 Color code: A, B
Dimensions: 7-7/8" x 3-3/4" x 15/16"



20Z202
**No 16
Hand Instruments
(Deep Well)**

\$9.50 Color code: A, B
Dimensions: 7-3/4" x 3-3/4" x 1"



20Z203
**No 17
Clamps, Retainers, etc.**

\$9.50 Color code: A, B
Dimensions: 7-7/8" x 3-3/4" x 3/4"



20Z204
**No 18
Finishing Strips**

\$9.50 Color code: A, B
Dimensions: 7-7/8" x 3-3/4" x 3/4"



20Z205
**No 19
Pliers, Scissors, etc.**

\$9.50 Color code: A, B
Dimensions: 7-7/8" x 3-3/4" x 3/4"



20Z206
**No 20
Discs, Wheels, Stones**

\$9.50 Color code: A, B
Dimensions: 7-7/8" x 3-3/4" x 3/4"



Lab Pans


Organizes and transports impressions, models, dentures, and crown and bridge appliances. Stainless steel spring clip holder. Color-coded.

17Z101
Lab Pan

\$7.75 Color code: A, B, D, G, M
Dimensions: 7-1/2" x 5-3/4" x 7"

17Z102
Wall-Hanging Lab Pan

\$7.75 Color code: A, B, D, G, M
Dimensions: 7-3/4" x 5-3/4" x 7"



17Z104
Wall Strip Lab Pan

\$10.50 White only
Drilled and countersunk with screws and wall anchors.
18" Wall Strip. Holds 3 Lab Pans.
Dimensions: 18" x 7/8" x 3/8"



E-Z Access Shelf & Accessories



Provides convenient work space for auxiliary items

20Z420 E-Z Access Shelf

\$130 Color code: A, G, I

E-Z mounting to standard equipment post (1-3/4" x 2-1/4"); No tools needed. Will not damage post's finish. E-Z adjustment to multiple positions Contoured edge to contain spills and materials. Stain resistant finish; Surface disinfect only. Weight capacity of 25lbs. Rotates 360°.

Standard E-Z Access Shelf as is accommodates 1-3/4" x 2-1/4" post.

Workable tray area: 13-3/4" L x 10-1/4" W



E-Z Access Shelf Adapters & Mounts S



20Z422 E-Z Access Shelf Adapter for Smaller Posts

\$7 Color code: A, G, I

Adapts to posts 1-3/8" to 2-3/4"

Dimensions: 2-1/8" x 1-7/8" x 3/4"

20Z437 E-Z Access Shelf 3 1/2" Bolt Adapter Kit

\$12 Color code: A, G, I

Accommodates post sizes from 2-1/4" to 2-3/8"

Dimensions: 4" x 2" x 2"

20Z438 E-Z Access Shelf 4 1/2" Bold Adapter Kit

\$12 Color code: A, G, I

Accommodates post sizes from 2-3/8" to 2-3/4"

Dimensions: 4-1/2" x 2" x 2"



20Z436I E-Z Access Shelf Cup Holder Snap-On

\$12 Gray only

Fits all standard cups. No tools needed.

Dimensions: 3-1/8" x 2-1/2" x 1"



20Z428 E-Z Access Shelf 8" Extension Arm

\$32 Color code: A, G, I

Enables shelf to extend over patient. Shelf has strength to hold 25lbs. when connected to arm. May be used with Wall Mount, Surface Mount, or directly to post.

Dimensions: 10" x 4" x 4-1/2"



20Z426 E-Z Access Shelf Wall or Cabinet Mount

\$32 Color code: A, G, I

Mount on vertical flat surface only.

Dimensions: 6-1/2" x 5-3/4" x 4-7/8"




20Z424 E-Z Access Shelf Surface Mount

\$17 Color code: A, G, I

Mount on horizontal flat surface only.

Dimensions: 3-3/4" x 3-3/4" x 2-7/8"



E-Z Storage S

E-Z Storage Tub System

Hinges allow cover to stay open on its own! Each tub comes with one vertical divider and one half-depth adjustment divider. Hinged cover comes in white or clear.

29R605
E-Z Storage Organizer (White Cover)

\$37 Rack sold separately

Dimensions: 5-7/8" x 4-7/8" x 5-3/4"

29R610
E-Z Storage Organizer (Clear Cover)

\$37 Rack sold separately

Dimensions: 5-7/8" x 4-7/8" x 5-3/4"

29R627
Replacement Hinges (Pair)

\$5.25

Dimensions: 7/8" x 1/2" x 1/8"

Hinges allow cover to stay open



Rack holding angled Tub

29R650
E-Z Storage Single Rack

\$23

Holds storage organizer Tub. Tub can be placed flat or at an angle for easy access of supplies. Wall or cabinet mounted with optional mounting bracket. Sturdy metal construction; plastic coated. White, non-slip feet.

Dimensions: 7" x 5" x 4-3/4"

29R651
E-Z Storage Mounting Bracket

\$5.25

White metal bracket for the E-Z Storage Tub single rack. Pre-drilled holes for cabinet or wall mount.

Dimensions: 5-1/8" x 1" x 1/2"



Handi-Hopper

Handi-Hopper (With Mounting Bracket) S

Small plastic waste receptacle for disposal of contaminated waste at chairside, keeping tray clean. Slotted mounting bracket (included) for easy removal and replacement. To use for infectious waste, order red (Biohazard) liners. To use for other waste, order clear liners.

20Z404
Handi-Hopper

\$22 Color code: G, I

Dimensions: 7" x 4" x 6-3/8"

20Z408
Handi-Hopper Bracket

\$4.50 Color code: G, I

Dimensions: 2-1/4" x 1/4" x 1"

Mount this mini waste receptacle within arm's-reach of chairside for convenient disposal of contaminated waste

Handi-Hopper Liners

\$25

20Z405 | 100/BOX
Red Biohazard

Single use

20Z410 | 100/BOX
Clear

Single use



22Z100
Waste-Drop Ring

\$15 Color code: A, G, I

Waste-Drop Ring for cabinet top. Fits 5" hole. Made of self-lubricating, easy-to-clean resin.

Dimensions: 5-5/8" Dia. x 3"

S



HVE & SE Accessories

Saliva Ejector Valve (with Shut-Off)

- Quick disconnect with swivel
- Rubber tip with protective screen
- Provides variable control from full on to full off



Saliva Ejector Valve - Lever



HVE Valve - Lever

23E365
Saliva Ejector Valve with Lever On/Off Control

\$35 Gray only

Dimensions: 3" x 5/8" x 3/4"

23E364
Saliva Ejector Valve Lever Only

\$10 Gray only

Dimensions: 5/8" x 3/16" x 5/8"

23E363
Saliva Ejector Replacement Tip

\$10 Gray only

Dimensions: 5/8" x 5/8" x 1"

23E300
HVE Valve (Lever Style)

\$45 Silver only

Dimensions: 3" x 5/8" x 11/16"



Saliva Ejector Screens

Single use

Dimensions: 3/8" Dia. x 3/8"

\$12

25PK | 50Z952S

\$82

250PK | 50Z952P

Needle Capper



Needle Capper

\$29

One-handed capping and recapping. Syringe is ready for use at convenient position in Needle Capper. Base of device is secured to tabletop or tray with adhesive strip (included). Needle is uncovered at the same time syringe is withdrawn from Capper. After use, syringe is replaced for recapping. Needle Capper available in Standard and Long sizes. Either size tube will fit standard base.

Standard Needle | 27R700

White only

Includes: base, one tube, adhesive tape.


Dimensions: 2-3/16" x 1-1/4" x 1-1/4"

Long Needle | 27R704

White only

Includes: base, one tube, adhesive tape.

Dimensions: 2-1/2" x 1-1/4" x 1-3/4"



Water Wise S

Water Wise vs. Traditional Water Bottle

Features	Water Wise	Water Bottle
Easy to clean	X	
Antimicrobial	X	
Easy-to-fill large mouth	X	
Quick release	X	
Level indicator coded for: Sterile and distilled water, Medicament	X	
Strong	X	X
Easy to grip	X	
Sealable lid	X	
No threads to strip	X	
Pressure relief valve	X	



23Z100
Water Wise Starter Kit
\$34

Holds 1 liter. Includes bottle, adapter, lid, neon yellow and neon blue bubbles (water level indicator), and tube.

23Z105
Water Wise Replacement Bottle
\$23

Holds 1 liter. Includes bottle, blue lid, and neon yellow and neon blue bubbles (water level indicator), and tube.

- 1 Quick release adapter Pressure relief valve
- 2 Water level indicator (bubbles)
- 3 Flat stable bottom with large opening for easy cleaning and filling

Dispensers S



28R800
Cotton Roll Dispenser
\$34 Color code: A, B

Heavy-duty, no-skid base. Dispenses up to three rolls at a time. Extra sturdy construction.

Dimensions: 4-3/16" x 2-1/16" x 4"



55Z200
Matrix Dispenser
\$28.50 Color code: A, B

One roll of 3/8" x 12 yards Matrix clear tape 0.002 gauge. Allows matrix to be contoured (u-shaped) as it dispenses. Stainless steel cutter allows for desired length. Disposable.

Dimensions: 2-1/8" x 2-1/4" x 5/8"



Pellet Dispenser
\$34 White only

Spring loaded – pushes pellets to the surface.

28R825 | For Small Pellets (sizes 3, 4, 5)
 Dimensions: 1-7/8" x 1-7/8" x 2-1/16"

28R850 | For Medium Pellets (sizes 0, 1, 2)
 Dimensions: 1-7/8" x 1-7/8" x 2-1/16"

Dispos-A-Screen Traps & Cuspidor Screens



Dispos-A-Screen Traps

Avoid contact with potentially hazardous material. Improves performance of suction system. Single use.

\$19 25pk -or- **\$82** 144pk

Mauve 1 ¼"

Adec, Cascade, and Performer.

51Z080S | 25pk

51Z080V | 144pk

Green 1 7/8"

Ampco, Belmont, Biotec, Chayes, Dabi Atlante, Dansereau, DCI, Den-Tal-EZ, Dentech, Forest, Knight, Pelton & Crane, Proma, and others.

51Z123S | 25pk

51Z123V | 144pk

Blue 2 1/8"

Adec models with black plastic canisters, plus newer models by Ampco, Belmont, DCI, Forest, Marus MDT, Schein.

51Z045S | 25pk

51Z045V | 144pk

Yellow 2 ¼"

Adec, model with metal canisters.

51Z040S | 25pk

51Z040V | 144pk

Coral 2 5/8"

Pelton & Crane Spirit Unit with Plastic Canister

51Z050S | 25pk

51Z050V | 144pk

Gray 2 ¾"

Midwest 210 Models.

51Z070S | 25pk

51Z070V | 144pk

Beige 2 7/8"

Den-Tal-EZ models with metal canister.

51Z060S | 25pk

51Z060V | 144pk

Cuspidor Screens

Single Use

\$19 25pk -or- **\$82** 144pk

Large Universal 2 ½"

Gray

Adec, Ampco, Belmont, Dabi Atlante, Dansereau, DCI, Den-Tal-EZ, Forest, Knight, Marus MDT, Proma, Ritter, Schein.

51Z205S | 25pk

51Z205V | 144pk


Small 1 5/8"

White

Adec, Cascade, and Performer cuspidor bowls.

51Z200S | 25pk

51Z200V | 144pk



Mighty Mixer Bowls, Spatulas, and Organizers

Mighty Mixer Spatulas

Self-lubricating for easy cleaning. Broad, medium flex blade. Ideal for alginates and plastic.



50Z501
Single Spatula
\$6.50 Each White only

Dimensions: 7-5/16" x 3/16" x 1-1/4"



50Z501ASTD
Assorted 4 Pack
\$22

One (1) of each color: J, N, P, and R

Dimensions: 7-5/16" x 3/16" x 1-1/4"



50Z523
Delrin Spatula 10pk
\$11

Flexible flat blade. Non-stick. Resistant to chemicals. Great for mixing cements, acrylics, etc.

Dimensions: 4-5/8" x 3/16" Dia.



50Z502
Cement Spatula
\$7 Each

White, plastic, double-end spatula. For fillings, cement, and pastes. Self-lubricating plastic for easy cleaning.

Dimensions: 6-3/8" x 1/4" x 3/8"



50Z521
Mixing Sticks 100pk
\$23 White only

Single use.

Dimensions: 4-1/2" long



20Z496V
Mixing Well 144pk
\$24 White only

Single use.

Dimensions: 4-1/2" long



Mighty Mixer Bowls 36pk
\$23

Disposable, eliminating cleanup time. Fits securely into mixing bowl without slippage. Saves wear on mixing bowl. Free-standing for use with or without mixing bowl. Capable for use with electric Alginate mixers. Single use.

50Z505 | White
50Z505ASTD | Assorted

Dimensions: 4-7/8" Dia. x 3-1/8"



20Z490B
Implant Organizer
\$28 Blue only

Quick and easy recognition. Individual compartments for implants. Shape of the organizer matches the arch of the mouth for implant placement. Separate compartments prevent implant appliances from mixing together. Lubricated plastic enables tooth number to be written in the designated area. Cover helps keep materials secure.

Dimensions: 6-7/8" x 4-5/8" x 5/8"



20Z495R
Restorative Organizer
\$33 Neon purple only

Compartments shaped like the arch of the mouth for veneers, inlays, composite, and onlays. Disposable Mixing Wells included for bonding agent and cement. Accessory area for brushes and materials. Lubricated plastic enables tooth number to be written in the designated area. Light protective cover can be used as a light shield and to prevent materials from setting early. Includes 25pk of Mixing Wells.

Dimensions: 6-7/8" x 4-5/8" x 5/8"



win the 😊 dental day

We know taking exceptional care of your patients is job number one, which is why every Zirc product is designed to help you reduce moments of high emotional stress and eliminate inefficiency within your routine tasks and procedures. Our clinically-trained consultants are ready to help you lead the way in transforming your practice into the most efficient version of itself. Learn more at zirc.com.

Who we are



Reliability



Passion



Creativity



Customer Focus





Integrity



Find Solutions

Made in
United States





ZIRC.COM | (763) 682-6636

**All U.S. pricing subject to change*

20Z500C REV B 0424