



CELEBRATING  
**50**  
YEARS

# Zirc

Because Time Is Everything

**PRODUCT GUIDE 2018**

[www.zirc.com](http://www.zirc.com) | 800.328.3899

## WHY ZIRC

2

## COLOR CODING

4 - 21

Why Color Coding? .....	4
Color Chart .....	5
Tub Kits .....	7
Tub Accessories .....	8
Trays & Covers .....	10
Multi-Rod Racks .....	12
Cassettes .....	13
Steri-Containers/Mats .....	15
Bur Guards/Blocks .....	16
E-Z ID Rings & Tape .....	19
Zobe .....	20

## ISOLATION

22 - 25

mr.thirsty® one-step .....	22
Insti-Dam® I SE Cushion I Lingua-Fix® .....	24

## MOUTH MIRRORS

26 - 31

Crystal HD® Mirrors .....	26
mirror magic® Anti-Fog Solution .....	30

## OPERATORY/MISC

32 - 40





## Zirc History 1967 - 2017

James Campion, our CEO, founded Zirc in 1967 and his daughter Nicolle, joined the team in 2000, and is the companies President. In 2017 we celebrated 50 years of offering our customers the highest quality products at the best value. Dental technology and processes have seen many changes over the last 50 years, but at Zirc, our commitment to providing time-saving, innovative dental products continues to be our top priority. We look forward to serving our customers for years to come!

## Zirc Perks

Members receive FREE samples, money-saving promotions, tips to get organized and stay organized, product evaluation invites, and more! Sign up to become a Zirc Perks Member via [www.zirc.com](http://www.zirc.com) or call us at +1 763-682-6636 (available in U.S.A. and Canada only).

## Made in the USA

Our headquarters, located in Buffalo, Minnesota, houses our engineering, tooling, manufacturing, distributing, and customer service operations. We're proud to say that our state-of-the-art, high-quality products are made in the U.S.A.

## Antimicrobial

Antimicrobial compounds are intended to protect parts against staining, odors, and degradation due to fungi, mold, yeasts, mildew, algae, and other microbes. Zirc incorporates antimicrobial compounds during the manufacturing process. Unlike a coating, this prevents fading or wearing away, giving you long-lasting protection. The incorporation of antimicrobial compounds in our plastic products is not intended to protect users or others against food-borne or disease-causing bacteria. For more information, please reference EPA Registration Notice 2000-1.



## Experience Our Products\* for FREE

The ZOBE kit contains over 30 of our top-selling products for a stress-free, fun, organized practice. ZOBE ships at no cost, allowing you to explore the products before you buy. No charge shipping and return. See page 20 for more details.

\*for demo/viewing only




## WANT TO MAXIMIZE SPACE?

**E-Z Jett® Cassette 8-place Slim**

Our autoclavable E-Z Jett® Cassettes extend the life of your instruments by protecting and keeping them together through your procedures and sterilization process. They are designed for superior water flow during cleaning and have non-slip feet for stability and safety.

- Holds 8 instruments
- Covers attached to the cassette open and fold back easily
- Compact size - Fits the Statim
- Maximizes the space of the Ultrasonic and Autoclave
- Stackable in Ultrasonic. Cassettes lightly attach to one another
- Antimicrobial incorporated during manufacturing process
- Syringe Tip Holder included
- Autoclavable

**Maximize space!  
Stackable in Ultrasonic  
and Autoclave!**



## FUN &amp; EASY PERSONALIZATION

**Imprinting Features Now Available**

We now offer new and improved personalized imprinting options for our popular denture and retainer boxes. Add your logo or unique message or choose one of our fun or professional designs.

**Choose Your Message with New Imprinting Options**

- Showcase Your Brand – Add contact info and upload your logo
- Make it Fun – From emojis to superheroes, our fun designs will bring a smile
- A Professional Image – Choose from one of our classic designs to enhance your brand


**Custom Imprinting!  
See page 32.**



## COMPACT &amp; CUSTOMIZABLE

**Performance Tool Organizer**

A deeper, roomier cassette made for specialty instruments that will not fit in standard cassettes. Design your own layout with our customizable, removable feet. The feet can be positioned to support and protect your valuable instruments. Use the Performance Tool Kit in conjunction with any of our E-Z Jett® Cassettes, which hold your hand instruments. Designed for superior water flow for superb cleaning. See page 15 for full details!







## ORGANIZE WITH EASE Why Color Coding?

### Increase Profitability

Color coding decreases time for setup and clean up, allowing you to see more patients. It also reduces the amount of inventory needed, limiting expired materials.

### Other Benefits

- Simplified training of new staff
- Decreases amount of inventory – fill one tub vs. stocking multiple treatment rooms
- Easy and efficient re-stocking of dental materials
- Reduces the amount of expired materials
- Reduces the risk of an exposure incident
- Eliminates the need to stop a procedure so the assistant can find a missing instrument or material
- Contains dirty instruments in covered, puncture-proof, leak-proof containers satisfying the CDC and OSHA standards
- Small items stay organized and at your fingertips
- Set-ups appear organized, sanitary, and visually pleasing to patients

**Want help on how to color code in your dental office?**

*For more information, visit us online to download our Color Coding brochure.*



**Color coding has been proven to save time for each procedure and will allow you to see more patients in a day.**



# ZIRC Color Chart

## Classic



White  
A



Blue  
B



Green  
D



Beige  
G



Gray  
I



## Jewel



Plum  
E



Teal  
J



Midnight Blue  
T



Gray  
I



Red  
M



Copper  
U



## NEW COLORS

## Vibrant



N. Blue  
N



N. Yellow  
O



N. Green  
P



N. Orange  
Q



N. Purple  
R



N. Pink  
S







Syringe Composite Tub



Capsule Composite Tub



Endo Tub



Crown Prep Tub



Hygiene Tub



Implant Tub

## SMALL CHANGES Big Results

**Cluttered drawers = stress and chaos. Stay efficient with color-coded organization.**

Start by removing everything from your operatory drawers. Place in a large box and set aside. Restock your operatory drawers with staple supplies such as: bibs, saliva ejectors, cotton products, etc. (basically, anything you use in EVERY procedure).

Here's where the Tubs come into play. You need a minimum of one Tub for each procedure. Each Tub should be a different color so that you can quickly identify procedures by color. Stock each Tub with all materials you will need for that procedure.

Store ALL Tubs in a central area (such as sterilization). Take the Tub into the treatment room when a procedure is performed. Pull out what you need at chairside. Disinfect and re-store.

**Reduce your stress and have fun  
with color organization.**





## MAKE IT SIMPLE WITH AN All-In-One Kit

Our kits are designed to store and organize materials and keep your inventory under control. In addition, they are easily portable for transportation between procedures and your centralized location for sterilization, disinfecting, and restocking.

### Complete Tub

Includes: Procedure Tub, Tub Liner, Tub Dividers, Cup/Cover Single, Cup/Cover Double, Divided Slide Tray, Safe-Lok Cover for Tub®. Dimensions: 12-5/8" x 11-1/4" x 4"

**\$66 (20Z455) Available in all colors**

AM S



Use: Any Procedure

## Composite Kits


Store and organize composite materials with our convenient Composite Kits. When stored in a Tub, contents become portable and can be transported to a central sterilization area for disinfecting and restocking. Organizational tools help you keep better inventory control with Clip-On Holders that allow you to separate composite for use during procedures. Included Long Tub Cup and Slide Tray provide additional organizational space and Shade Label Sheet provides self-stick labels preprinted with composite shade.

### Syringe Composite Kit

Syringe holder allows you to store up to 40 syringe tubes. Kit includes: Syringe 20-Unit Stand, 1 Long Tub Cup/Cover, 1 Divided Slide Tray, 1 Syringe Clip-On Holder, Shade Label Sheet 2pk.

**\$52 (20Z489)**

AM S



### Capsule Composite Kit

Long Tub Cups provide a total of 20 compartments for capsule storage. Includes: Tub Liner, 4 Long Tub Cup/Cover, 1 Divided Slide Tray, 1 Capsule Clip-On Holder, Shade Label Sheet 2pk.

**\$52 (20Z479)**

AM S





## Build Your Own Tub



### Procedure Tub

Allows you to “build” your own setup.  
Dimensions: 12-1/4" x 10-7/8" x 2-3/4"

**\$19** (20Z463) Available in all colors



### Safe-Lok® Tub Cover

Locking cover for securing materials.  
Keep materials protected during storage and transportation.  
Dimensions: 12-1/4" x 10-7/8" x 2-3/4"

**\$18.50** (20Z459)



### Tub Liner

Foam liner that sits in bottom of Tub.  
Tub not included. Replace as needed.  
Dimensions: 9-5/8" x 8-5/8" x 1/8"

**\$6** (20Z462)



### Tub Box

Designed to hold/organize “hard-to-find” supplies.  
Non-slip bottom for secure placement.  
Dimensions: 3-3/4" x 3-5/8" x 1-7/8"


**\$15** (20Z474) Gray Only



### Tub Insert

10 compartments for easy organization.  
Sits on rim of Tub, allowing for storage underneath.  
Dimensions: 10-3/4" x 9-11/16" x 1"

**\$14.50** (20Z465) Gray Only



### Tub Dividers

Creates up to 24 compartments.  
May be cut or trimmed to create larger compartments.  
Dimensions: 10-1/16" x 9-1/16" x 1-5/16"

**\$12** (20Z461) Gray Only



Use color-coded tape on Tub Cups and Tub Box to identify procedural contents. (Learn more on page 19)



### Tub Cups with Covers

Fill Tub Cups with brushes, clamps, cotton pellets, prophylaxis angle tips, prophylaxis paste, pumice, syringe tips, implant tools, etc. Can be used in procedure tubs, drawers, on trays, etc. for small, easy-to-lose items. Covers are transparent.



### Single Tub Cup with Cover

Dimensions: 2-1/4" x 1-5/8" x 1-1/8"  
**\$7.25** (20Z471) Gray Only

### Double Tub Cup with Cover

Dimensions: 3-3/8" x 2-1/4" x 1-1/8"  
**\$7.25** (20Z472) Gray Only

### Long Tub Cup with Cover

Dimensions: 8-3/8" x 2-5/16" x 1-1/8"  
**\$13.25** (20Z473) Gray Only

# Tub Accessories



## Divided Slide Tray

Dimensions: 9-3/4" x 5-3/8" x 1-1/16"

**\$11.50** (20Z466) Clear

## Flat Slide Tray

Dimensions: 9-3/4" x 5-3/8" x 1-1/16"

**\$11.50** (20Z467) Clear



## Shade Label Sheets 2pk

Preprinted labels for marking Tub materials by procedure. Shade Labels can be used with the Syringe Composite and Capsule Composite Tub Organizers.

Dimensions: 7" x 4-3/8" x 1/64"

**\$5** (20Z477)

## Slide Trays

Sits on rim of Tub and allows for storage underneath.



*Check out our Implant and Restorative Organizers on page 40.*



## Syringe 20-Unit Stand

Stands up for easy access to composite.

Lays flat to fit inside Tub or drawer.

Holds up to 20 syringe tubes.

Has space for 2 Syringe 10-Unit Add-Ons (Stand can hold up to 40 syringes).

Dimensions: 8-3/4" x 4-3/4" x 1-9/16"

**\$25** (20Z487)



## Syringe 10-Unit Add-on

An Add-on for the Syringe 20-Unit Stand.

Gives additional space for syringe storage.

Holds 10 syringe tubes.

Dimensions: 8-9/16" x 2-5/16" x 3/4"

**\$10** (20Z488)



Syringe 20-Unit Stand

*The Syringe 10-Unit Add-on may be added for storage of up to 40 syringe tubes.*



## Syringe Clip-On Holder

Allows quick access to composite during procedure.

Conveniently clips on any side of Tub, B-Lok Tray or Capsule Composite.

Holds 15 syringe tubes and 3 applicator tips.

Dimensions: 8" x 2-7/16" x 3/4"

**\$10** (20Z486)



## Capsule Clip-On Holder

Allows quick access to composite during procedure.

Conveniently clips on any side of Tub or B-Lok Tray.

Holds 8 capsules, 4 syringe tubes, and capsule/gun dispenser.


Dimensions: 8" x 2" x 9/16"

**\$10** (20Z476)





## Trays & Covers



### Are you in compliance?

#### Trays

Our exclusive Safe-Lok® tray cover provides an excellent method of infection control and safety that allows for two benefits other trays do not:

- 1 Prevents any contaminated items from falling off of the used tray during transportation from operator to central sterilization. The contaminated items stay enclosed until you are ready for the cleaning and sterilization process.
- 2 Set up your procedure trays ahead of time with all instrumentation needed. The Safe-Lok® cover securely locks to the tray, not only keeping items clean, but also allowing for safe transportation to chairside. The cover is removed once the patient is in the chair and the procedure begins.

#### OSHA

Immediately or as soon as possible after use, contaminated reusable sharps shall be placed in puncture-resistant, leak-proof and labeled or color-coded containers until properly reprocessed. (To avoid spillage of contents, it is suggested that they be covered and secured prior to moving.)

#### CDC

Dirty Instruments must be "contained" prior to transporting them. Containment can be defined as **instrument cassettes** or cages in which dental care items are cleaned, sterilized, and stored until point of use. Sharp instruments should not be carried openly to the instrument-processing area. Percutaneous injuries may occur to other dental healthcare personnel (DHCP) or patients in hallways leading to the instrument-processing area.

**The Safe-Lok® Tray Cover helps prevent cross-contamination and provide safe transportation back to the sterilization area.**

# Trays & Covers



## B-Lok Flat Tray

Wide edges for secure grip.  
Non-slip feet for stability.  
Self-lubricating, high-gloss surface helps prevent staining.  
Dimensions: 13-3/8" x 9-5/8" x 7/8"

**\$12.50** (20Z401) Available in all colors



## B-Size Safe-Lok® Tray Cover

Fits B-Lok Flat and Divided Trays only.  
Side locking mechanism securely locks cover to tray.  
Dimensions: 13-7/8" x 9-7/8" x 3/4"

**\$17.50** (20Z445)




## B-Lok Divided Tray

Holds 12 hand instruments.  
Bur holder for 9 FG, RA, or HP burs.  
Comes with Insti-Grip.  
Non-slip feet for stability.  
Self-lubricating, high-gloss surface helps prevent staining.  
Holder for Syringe Tip or Saliva Ejector.  
Dimensions: 13-3/8" x 9-5/8" x 7/8"

**\$14.50** (20Z451) Available in all colors



## FOR B-LOK DIVIDED TRAY



### Insti-Grip

To be used with the B-Lok Divided Tray only.  
Holds instruments quietly and safely in place.  
Dimensions: 6-5/8" x 1/4" x 3/8"

**\$2** (20Z454I)



## Mini Tray

Self-lubricating, high-gloss surface helps prevent staining.  
Non-slip feet for stability.  
Wide edges for secure grip.  
Mini Trays are great for hygiene/exam setup.  
Dimensions: 9-3/8" x 6-3/8" x 7/8"

**\$11.50** (20Z101) Available in all colors



## Mini-Lok Tray Cover

Side locking mechanism securely locks cover to tray.  
Dimensions: 9-7/8" x 6-5/8" x 1"

**\$15** (20Z103)



*Our E-Z Jett® Cassettes fit on Zirc Trays with locking cover! See page 13 for E-Z Jett® Cassettes.*



**The most efficient method  
of infection control.**

## B-Size Tray Cover (Non-Locking)

Dimensions: 14" x 10-1/8" x 5/16"

**\$17** (20Z441)



## Mini Tray Cover (Non-Locking)

Dimensions: 9-7/8" x 6-5/8" x 1"

**\$14.50** (20Z102)





# Multi-Mod Racks

## SAVE ON SPACE & Organize Your Office

Tub, tray, and cassette racks are very effective additions to a dental office. They save on counter space and several tubs or trays with covers can fit into one rack. These racks have non-slip feet so they are very stable on any counter or cabinet surface. With a chrome-plated finish, these racks can fit any dental office design.

### Multi-Mod Rack 8-Place

Holds 8 Trays or 4 Tubs with covers.  
Adjustable to hold Trays, Tubs, and Cassettes.  
Sturdy metal construction; chrome-plated.  
Non-slip feet.

**\$92** (21Z410)

S

### Multi-Mod Rack 6-Place

Holds 6 Trays or 3 Tubs with covers.  
Adjustable to hold Trays, Tubs, and Cassettes.  
Sturdy metal construction; chrome-plated.  
Non-slip feet.


**\$92** (21Z105)

S

Multi-Mod Rack (8-Place)



Multi-Mod Rack (6-Place)



**Adjustable  
to fit your needs!**



*Can adjust to hold  
Mini Trays and Cassettes*





## WHY USE Zirc's Cassettes?

A low-cost investment with a quick return! Individual instrument handling takes time away from treating patients and increases overall processing time. Zirc's cassettes will prolong your instrument's life and save you time without breaking the bank.

In a 2014 Implementation of Cassettes Case Report, Marie T. Fluent, DDS, writes how the dental office benefited from the transition from no cassettes to using them with every procedure: "Drs. Dale and Thomas were both grateful that they chose to implement ... an Instrument Management System for the organization, processing, and handling of their instruments within their practice. They agreed that cassettes improved all aspects of instrument management... Initially, their opinions were varied; however, at the conclusion of the one-year anniversary, no clinical staff could imagine processing instruments in any other fashion. Cassette usage had become organized, efficient, safe, and routine." Read the article in its entirety at: [www.dentaleconomics.com](http://www.dentaleconomics.com).

- Designed for superior water flow and cleaning - open top, round holes, and radius corners
- Keeps instrument setups together through sterilization and procedure processes
- Protects instruments from extensive damage and prolongs their life
- Non-slip feet for stability and safety
- AUTOCLAVABLE
- Available in 16 colors for color coding by procedure

Keep the workplace  
organized and confusion  
to a minimum!



Learn more at  
[www.zirc.com/instrumentcassettes](http://www.zirc.com/instrumentcassettes)

## Sterilization Pouch Size

Cassette	Pouch Size
E-Z Jett® 5-Place	5-1/4" x 12"
E-Z Jett® 8-Place	7-1/8" x 13-13/16"
E-Z Jett® 8-Place Slim	7-1/8" x 13-13/16"
E-Z Jett® 10-Place	10" x 15"

*Provide pleasing presentation chairside.  
Minimize stress for patients.*

## Sizing Chart for E-Z Jett® Cassette Sterilization

Autoclave Manufacturer	Model	Chamber Dimensions L x W x D	E-Z Jett® 5-Place	E-Z Jett® 8-Place	E-Z Jett® 8-Place Slim	E-Z Jett® 10-Place
Midmark	M11	11" x 18"	12	8	8	4
	M9	9" x 15"	6	4	4	2
Pelton & Crane	Magnaclave	15-1/4" x 15" x 26"	26	24	24	12
	Delta XL 10	9-13/16" x 17-1/2"	12	9	9	4
	DeltaXL8	8-1/4" x 14"	8	4	4	2
SciCan	Bravo17	10" x 13-1/2"	9	6	6	3
	Bravo21	10" x 17-1/2"	12	9	9	6
Tuttanuer	3870	15" x 27"	20	12	12	10
	EZ10	10" x 18"	12	9	9	6
	2540M	10" x 18"	12	9	9	6
	EZ9	9" x 18"	8	6	6	4

## AUTOCLAVABLE Instrument Holders




Syringe Tip Holder

**Syringe Tip Holder 5pk**  
Dimensions: 1" x 3/8" x 3/8"  
**\$5 (50Z926)**

### E-Z Jett® Cassette 5-place

Holds 5 instruments.  
Great for a basic exam set up!  
Syringe Tip Holder included.  
Covers stay attached to the cassette and fold back easily.  
Fits Statim Sterilizer.  
Dimensions: 8" x 3-1/8" x 1-3/16"

**\$27 (50Z925) Available in all colors**



### E-Z Jett® Cassette 8-place

Holds 8 instruments.  
Covers stay attached to the cassette and fold back easily.  
Syringe Tip Holder included.  
Dimensions: 7-15/16" x 4-7/16" x 1-1/8"

**\$37 (50Z927) Available in all colors**



### E-Z Jett® Cassette 8-place Slim



Holds 8 instruments.  
Covers stay attached to the cassette and fold back easily.  
Compact size. Fits the Statim.  
Stackable in Ultrasonic.  
Syringe Tip Holder included.  
Dimensions: 7-1/16" x 4-1/16" x 1-5/16"

**\$32 (50Z911) Available in all colors**



### NEW PRODUCT


Size comparison with the  
8-Place E-Z Jett® Cassette.  
See page 3 for more details.



### Instrument Mat

Holds 10 hand instruments.  
Reverse side holds 4 surgical instruments.  
Non-slip, soft material.  
Small size allows easy access to instruments.  
Organizes instruments in order of procedure use.  
Can be cut to fit Mini-Tray.  
Dimensions: 6-3/4" x 2" x 3/8"

**\$14 (20Z464) Available in all colors**



Instrument Strap

### Mini Mat


Holds 8 instruments.  
Organizes instruments at chairside and with sterilization.  
Non-slip feet for stability.  
Dimensions: 3-7/8" x 3-1/2" x 1/2"

**\$12.50 (20Z480) Available in all colors**

### Instrument Strap 2pk

Locking strap secures instruments.  
Dimensions: 3-5/8" x 1/2" x 5/16"

**\$6/pair (50Z917) Gray Only**



### E-Z Jett® Cassette 10-place

Holds 10 instruments.  
Accessory compartment in the middle for extra storage.  
Covers stay attached to the cassette and fold back easily.  
Dimensions: 7-15/16" x 7" x 1-1/8"

**\$42 (50Z930) Available in all colors**



# Steri-Containers



Closed

## Performance Tool Organizer

A deeper, roomier cassette made for large and specialty instruments. Design your own layout with our customizable, removable feet (also sold separately). Designed for superior water flow for superb cleaning. Compact styling to fit in most sterilizers. Covers stay attached to the cassette and fold back easily. Ideal for:

- Surgical elevators and forceps
- Orthodontic pliers
- Hemostats
- Scissors
- Any instrument that does not fit in your regular cassette



Open

### Performance Tool Organizer

Covers stay attached to the cassette and fold back easily.  
Dimensions: 8-3/16" x 4-1/4" x 1-7/32"

**\$55 (50Z940) Available in all colors**



### Performance Tool Organizer Replacement Feet

Two strips of 10.

**\$7 (50Z944) Gray Only**



**Holds your large and specialty tools**

*E-Z Jett® Cassettes fit on Zirc Trays with locking cover! Check out our Trays & Safe-Lok® Cover on page 10.*



## Steri-Container

Holds 10-14 instruments.  
Keep instrument setups together through sterilization processes.  
Dimensions: 8-1/8" x 1-7/8" x 1-7/8"

**\$11.50 (50Z900) Available in all colors**



## Compact Steri-Container

Holds 8-12 instruments.  
Easy open latch, just push down.  
Statim 2000 holds 7 Compact Steri-Containers.  
Dimensions: 7-1/8" x 1-1/2" x 1-1/2"

**\$11.50 (50Z905) Available in all colors**





## Bur Guards



### PROTECT YOUR BURS Steri-Bur Guards




### Steri-Bur Guards and Adapters

This gadget has it all! Not only does it protect your burs during cleaning and sterilization, but it also can be positioned at three different angles to give you easy access at chairside.

#### Steri-Bur Guard

- Fits FG, RA, and HP Burs
- Adjustable height for multiple-length burs (12-Hole Tall Bur Guard not adjustable)
- Holds bur in position for maximum cleaning effectiveness
- Drain holes help prevent rusting of burs



#### 12-Hole Bur Guard

Accommodates burs from 12mm-30mm.  
Dimensions: 2-7/8" x 3/8" x 1-3/8"

**\$15 (50Z406) Available in all colors**



#### 22-Hole Bur Guard

Accommodates burs from 12mm-30mm.  
Dimensions: 4-7/8" x 7/16" x 1-7/16"

**\$16 (50Z410) Available in all colors**



#### 12-Hole Tall Bur Guard

Accommodates burs from 30.9mm-37.1mm.  
Dimensions: 2-7/8" x 3/8" x 1-5/8"  
Ideal for nose-cone lab burs


**\$15.50 (50Z408) Gray Only (I)**

*A Short Bur Adapter is included to stop any shorter burs from falling out.*

#### Universal Bur Adapter 5pk

Can be used with 12- or 22-Hole Bur Guard.  
Enables long and short burs to be used in the same bur guard without short burs falling out.

**\$5 (50Z412)**



# Endo Guards & Holders




## Endo Guards

A unique delivery system to organize, store, and sterilize endodontic files and reamers. The Endo Guard swivels 360° with stops at every 45° for operator's choice. Large slotting provides complete visibility and draining during cleaning.

### Steri-Endo Guard


- Holds most lengths for staging, cleaning, and sterilizing files and reamers
- Millimeter gauge measures and positions rubber stopper
- Counter on locking bar tracks the number of sterilizations
- Right- or left-hand use; optional angling
- Available in Hand or Engine styles
- Single-use foam insert cleans/holds files (packaged with 12)



### Steri-Endo Guard Engine

Organizes up to 11 engine and 5 hand files/reamers.  
Dimensions: 5-3/8" x 3/8" x 2-3/16"

**\$30 (50Z452) Available in all colors**



### Steri-Endo Guard Hand

Organizes up to 16 hand files/reamers.  
Dimensions: 5-3/8" x 3/8" x 2-3/16"

**\$30 (50Z450) Available in all colors**

### Foam Inserts

Single use.

**\$30 (50Z451) 48pk**

**\$60 (50Z451V) 144pk**



*Unique organizer rests on your finger.*

### Endo Assist

Built-in millimeter gauge measures and positions rubber stopper. For right- or left-hand use. Small and large rings are included to accommodate most finger sizes. Flat bottom enables the Endo Assist to stand freely on its own. Single-use foam insert cleans/holds files (packaged with 12). Dimensions: 1-1/2" x 1-3/4" x 1-5/8"


**\$28 (50Z460) Available in all colors**



### Assist Stand

Dense foam cleans bio-debris faster and more easily. Holds Endo Instruments in an upright position. Square shape fits great in procedure tub. Thicker foam prevents solution from leaking. Rounded corners for easy cleaning. Dimensions: 1-7/8" x 1-7/8" x 1-3/8"

**\$16.50 (50Z462) Black Only**



### Assist Stand Foam

Single use.

Dimensions: 10-1/2" x 7-3/4" x 7/8"

**\$31 (50Z466) 48pk**

Dimensions: 10-1/2" x 7-3/4" x 7/8"

**\$61 (50Z466V) 144pk**

## Bur Blocks



### Magnetic Bur Blocks

- Fits FG, RA, and HP Burs
- Magnetic base holds burs securely
- Non-slip feet for stability
- Drain holes help prevent rusting of burs



#### 8-Hole Bur Block

Covers sold separately.  
Dimensions: 1-1/2" x 1-1/2" x 1/2"

**\$14 (50Z400)** Available in all colors



#### 14-Hole Bur Block

Covers sold separately.  
Dimensions: 2-15/16" x 1-1/2" x 1/2"

**\$15 (50Z403)** Available in all colors

Recessed area to hold used burs!



S

### Bur Block Covers

- Covers available to maintain cleanliness of burs during storage
- 1" covers for FG and RA burs; 2" cover (14-hole only) for HP lab burs
- Covers are not autoclavable



#### 8-Hole Bur Block Covers 1" High - Standard

Dimensions: 1-5/8" x 1-5/8" x 1-1/16"  
**\$7 (50Z401)**



#### 14-Hole Bur Block Covers 1" High - Standard

Dimensions: 3-1/8" x 1-5/8" x 1-1/16"  
**\$7 (50Z404)**



#### 14-Hole Bur Block Covers 2" High - HP Lab Burs

Dimensions: 3-1/8" x 1-5/8" x 2-1/4"  
**\$7 (50Z405)**



# Rings & Tape



## How to effectively identify your instrument setup.

### 1 color for practitioner:

Neon Orange shows ownership to the practitioner.

### 1 color for procedure:

Neon Blue is for identifying the procedure.

## Did you know?

*The diagonal setup of the orange Rings shows the order of how the instruments should be used. This eliminates confusion between the doctor and assistant.*

Place your instrument to determine which ring size you need.




Now available in XL size!



## E-Z ID Rings

- Prepackaged in an easy-to-use dispenser
- Available in Small, Large, and XL sizes to fit all instrument handles



### 25pk

Available in all colors.

- SM (70Z100) \$7
- LG (70Z200) \$7.25
- XL (70Z000) \$7.50

### Jewel System 200pk

Assorted colors (E, J, T, I, M, and U).

- SM (70Z115) \$49
- LG (70Z215) \$50
- XL (70Z015) \$51

### Vibrant System 200pk

Assorted colors (N, O, P, Q, R, and S).

- SM (70Z110) \$49
- LG (70Z210) \$50
- XL (70Z010) \$51

### Classic System 200pk

Assorted colors (A, B, D, G, and I).

- SM (70Z105) \$49
- LG (70Z205) \$50
- XL (70Z005) \$51



## E-Z ID Tape

- Tape dispenser has built-in tape cutter.
- Dispenser roller prevents tape from twisting.



### E-Z ID Tape Roll

Available in all colors.  
Includes : 1 - 10' roll.

\$10 (70Z300)



### E-Z ID Tape System - Jewel

Assorted colors (E, J, T, I, M, and U).  
Includes : 8 - 3' rolls.

\$42 (70Z315)

### E-Z ID Tape System - Vibrant

Assorted colors (N, O, P, Q, R, and S).  
Includes : 8 - 3' rolls.

\$42 (70Z310)

### E-Z ID Tape System - Classic

Assorted colors (A, B, D, G, and I).  
Includes : 8 - 3' rolls.

\$42 (70Z305)

Label Containers!  
Keep containers organized  
by labeling them per content.



Zirc's **O**rganizational **B**ox of **E**fficiency

### Why Z.O.B.E.?

- The ZOBE kit ships to your practice at NO cost – Free ZOBE, free shipping, free return
- ZOBE contains over 30 of our top-selling products for organization
- Have fun exploring the products – **Please Note: Z.O.B.E. is for demonstration purposes ONLY**
- Your ZIRC efficiency expert consults with you and listens to your needs
- Your ZIRC expert reviews the amount of patient flow for each procedure
- You receive a custom plan that's right for you
- Virtual consultations available

Visit us at [www.zirc.com/zobe](http://www.zirc.com/zobe)

Dreaming of a stress-free,  
fun, and organized practice?  
**Let us send you Z.O.B.E.**

\*Available in USA only.

**LEARN MORE 800.328.3899**

The ZOBE kit ships to your practice at NO cost.\*  
Free ZOBE. Free Shipping. Free Return.



**NEW!**  
Organize your  
practice with  
RESET!

# RESET

We come to you. We get you organized.

## HIT RESET AND TAKE CONTROL

Zirc and the Academy of Chairside Assisting (ACA) have partnered together offering to RESET your practice and get a fresh start with dental efficiency. Having one of our EFFICIENCY CONSULTANTS come into your practice can reduce stress, increase productivity and meet infection control guidelines.

Contact us to learn more: 800.328.3899





Pedo/Small Adult

Medium/Large



## The Efficiency of Hands-Free Isolation

A hands-free isolation device that retracts, isolates, and evacuates! Work continuously without the repeated use of suction, minimize your isolation inventory, and gain a quick return. With two convenient color-coded sizes, these one-time use isolation devices are the perfect cost-effective tool for your dental office.

- Built-in bite block for patient comfort
- Connects directly to HVE Valve
- Keeps work environment dry

### Pre-Assembled

The new mr. thirsty® one-step is ready to use and is single use! Just insert into your HVE Valve (one-step!).

### 2 Sizes Available

mr. thirsty® one-step comes in two convenient sizes, Pedo/Small Adult and Medium/Large. It can also be trimmed for a custom fit!



### 10pk Trial Kit

Includes 5 of each size:

- (5) Pedo/Small Adult
- (5) Medium/Large

**\$35** (50Z987)



### 25pk

- Pedo/Small Adult **\$70** (50Z986SM)
- Medium/Large **\$70** (50Z986LG)

ALSO  
AVAILABLE IN  
100 PACK



Replacement HVE  
Valve can be found  
on page 33

### 100pk

- Pedo/Small Adult **\$225** (50Z986SM-B)
- Medium/Large **\$225** (50Z986LG-B)

### Comfort Kit

The Comfort Kit's lightweight tubing allows your heavy HVE Valve/Tubing to stay mounted to the equipment.

Includes:  
Lightweight lever HVE Valve  
6' lightweight tubing  
mr. thirsty® one-step NOT included  
**\$60** (50Z992)



Take comfort to the next level with lightweight tubing!

mr. thirsty® one-step benefits everyone! Here's how...



**Doctor**  
**25% average decrease** in procedure time. Hooks directly to HVE Valve. Single-Use.



**Hygienist**  
**Reduce chair time by as much as 20%.** Hands-free continuous suction. Ease patient anxiety.



**Assistant**  
**Free up time.** See another patient in the next room. Catch up with sterilization.



**Patient**  
**Improved patient comfort.** Non-threatening appearance and improved patient safety.




**REDUCE CHAIR TIME BY AS MUCH AS 20%**



How does mr. thirsty® one-step stack up against the competitors?

mr. thirsty® one-step Competitor Brand

Cost-Effective	✓	
Large Upfront Investment		✓
No Maintenance or Sterilization Required	✓	
Can be Trimmed to Fit	✓	
Colored Product for Oral Cavity Contrast	✓	





To read a full product evaluation visit...  
<http://www.dentalproductshopper.com/relaxed-fit-insti-dam>


**Easy to use, non-threatening,  
 and comfortable for patient.**



## Insti-Dam®

### Features for BOTH Relaxed Fit and Original Insti-Dam®

- Built-in flexible medium gauge material that comes pre-assembled with frame and pre-punched hole
- Pre-punched hole makes positioning easy and allows eyes and nose to be open and comfortable
- Helps eliminate tearing (additional holes can be punched)
- Compact design fits outside patients lips, easy to use, non-threatening and comfortable for patient
- Radiographs may be taken by folding it to the side
- Single use



### Latex-Free Insti-Dam® Relaxed Fit 20pk

Ideal for posterior isolation.  
 The extra "relaxed" material in the Relaxed Fit Insti-Dam® makes for quick posterior placement - no stretching, tearing, or pulling.  
 Dimensions: 4-7/8" x 4-7/8" x 6-1/8" diameter

**Patent #7044738**

**\$70 (50Z457)**

### Blue Latex-Free Insti-Dam® 20pk

Ideal for anterior isolation.  
 Made with polyisoprene.  
 Dimensions: 4-1/4" diameter

**\$60 (50Z459)**

### Natural Latex Insti-Dam® 20pk

Ideal for anterior isolation.  
 Dimensions: 4-1/4" diameter

**\$55 (50Z455)**





## Lingua-Fix® & SE Cushion

### Lingua-Fix® Saliva Ejector

Lingua-Fix® provides a dry work area.  
Holds tongue in a steady position to protect from injuries.  
Simple operation... easy to place and remove.  
Comfortable, no sharp edges or corners.  
Large suction surface avoids blockage and is more comfortable.  
Reduces sedimentation of drilling by-products in pipes and tubing.  
Single use only.  
Comes pre-assembled.

**\$13.50 (50Z953S) 10pk**

**\$49 (50Z953) 50pk**



### SE Cushion

Used in conjunction with the Saliva Ejector.  
Prevents tissue grabbing and tissue irritation.  
Increases patient comfort.  
Improves range of suction.  
Quiets the noise of the suction.  
Single use.

**\$22 (50Z970) 100pk**



**HELPFUL TIP:**  
Ideal for use  
with rubber  
dams

**For comfortable and efficient  
suction or evacuation.**

### The Relax Factor:

The extra "relaxed" material in the Relaxed Fit Insti-Dam® makes for quick posterior placement - no stretching, tearing or pulling.



### Insti-Dam® Dispenser

#### Made for Insti-Dam® Original ONLY

Can be mounted to wall or cabinet or can sit on a counter.  
Available in White (A) for latex or Neon Blue (N) for latex-free.  
Non-slip rubber bottom for stability.  
Holds 35 Insti-Dams®.  
Dimensions: 4-7/8" x 4-7/8" x 6-1/8"

**\$20 (50Z471) Color Code: A & N**

AM S



A TV without HD is like a...  
A MIRROR WITHOUT HD



## Crystal HD


Experience high-definition quality. The brighter the workspace, the better the illumination and visibility. With ZIRC's Crystal HD® lens, you are visually working in a 40% brighter environment compared to the Rhodium lens that is currently the industry standard. If you haven't made the change to HD - what are you waiting for?

**40%  
BRIGHTER**  
Than Rhodium  
Mirrors



Your mirror **VS** Crystal HD®

### CUTTING EDGE IMAGERY



Rhodium = Reflects color 75% accurately  
Crystal HD® = Reflects color 99% accurately

Provides exceptionally bright,  
shadow-free images!





Crystal HD®

Available in multiple styles!

Designed to reduce eye strain and fatigue.

Autoclavable

”

*As a microscope user, I can understand the importance of clear vision. ZIRC mirrors significantly improved the brightness of my image, giving my eyes highest quality vision. It's true – brightness makes all the difference.*

- Dr. Dimitris Tsanaksidis, Microscopic & Aesthetic Dentistry

## CRYSTAL HD® MIRROR + MIRROR MAGIC® The Perfect Pair

Pair our Crystal HD® lens with our mirror magic® System for the highest quality vision and best procedure outcome. See page 30.

Crystal **HD** +   
mirror magic®  
anti-fog system





# Crystal HD® Mouth Mirrors & Handles

Resin

Stainless Steel

	RESIN	STAINLESS
Autoclavable	X	X
Antimicrobial protection	X	
Lens flush with mirror frame	X	
No galvanic shock on patient	X	
Rust-free	X	

**Why use Zirc's Resin Mirrors?**



## Assorted Soft Grip 12pk

### Vibrant Assorted

2 of each color: N, O, P, Q, R, and S.

**\$76 - #4** (50Z363-NEON)

**\$76 - #5** (50Z364-NEON)

### Jewel Assorted

2 of each color: E, J, T, I, M, and U.

**\$76 - #4** (50Z363-JWL)

**\$76 - #5** (50Z364-JWL)



## Assorted Thin Grip 12pk

### Vibrant Assorted

2 of each color: N, O, P, Q, R, and S.

**\$76 - #4** (50Z353-NEON)

**\$76 - #5** (50Z354-NEON)

### Jewel Assorted

2 of each color: E, J, T, I, M, and U.

**\$76 - #4** (50Z353-JWL)

**\$76 - #5** (50Z354-JWL)



## Black/Gray Soft Grip 12pk

**\$76 - #4** (50Z363)

**\$76 - #5** (50Z364)

## Black/Gray Soft Grip 4pk

**\$28.50 - #4** (50Z361)

**\$28.50 - #5** (50Z362)



## Handles



## Soft Grip Handle

**\$16 ea** (50Z340) Available in all colors



## Ergo Grip Handle

**\$16 ea** (50Z342) Color Code: B-S, and K (black)





## Resin Screw-on Heads 12pk

Flush surface, no debris build-up can occur.

### 1-Color 12pk

Available in all colors and K (black).

**\$76 - #4 (50Z371)**

**\$76 - #5 (50Z372)**

### Vibrant Assorted

2 of each color: N, O, P, Q, R, and S.

**\$76 - #4 (50Z371-NEON)**

**\$76 - #5 (50Z372-NEON)**

### Jewel Assorted

2 of each color: E, J, T, I, M, and U.

**\$76 - #4 (50Z371-JWL)**

**\$76 - #5 (50Z372-JWL)**

### Classic Assorted

A, B, D, G, I, and K (black).

**\$76 - #4 (50Z371-PSTL)**

**\$76 - #5 (50Z372-PSTL)**



## Vibrant Assorted Double-Sided Mirror Head 4pk

1 of each color: N, O, R, and S.

**\$40 - #4 (50Z374-NEON)**



## Stainless Steel Head 12pk

**\$76 - #3 (50R379)**

**\$76 - #4 (50R380)**

**\$76 - #5 (50R381)**



## Gray Double-Sided Mirror Head 4pk

**\$40 - #4 (50Z374)**



## Stainless Steel Head #0

**\$12.50 ea - #0 (50R378)**



**ZIRC's Mirrors & Handles are autoclavable!**

## mirror magic® Anti-Fog Solution



without mirror magic®



with mirror magic®



**Anti-Fog  
Solution removes  
debris from dental  
instruments!**

### Increase Your Chairside Efficiency

The mirror magic® Anti-Fog System prevents fogging on mouth mirrors and offers the practitioner a quick method for removing debris during procedures. The disposable Swipe Pads stick to the back of a gloved hand or a patient bib for easy access. The Anti-Fog System requires only a small amount of solution and is economic in size, making it easy to store.



Watch the magic in action  
at [www.zirc.com/mirrormagic](http://www.zirc.com/mirrormagic)



### mirror magic® Trial Kit

2-6ml bottles of Anti-Fog Solution  
24 Swipe Pads

**\$25 (50R200)**

Combine mirror magic® with Crystal HD®  
Mouth Mirror for the best view in dentistry.







*"The mirror magic® Anti-Fog System decreases fogging while also decreasing the build-up of debris on my mirror. This affords me outstanding visualization, which is crucial to performing quality dentistry. mirror magic® has truly made a difference in my practice and is one of those features I didn't realize I needed; and now I can't imagine practicing without it."*

— Andrew M. Goldsmith, DDS, DICOI, FIALD



## 10pk Anti-Fog Solution (6ml ea.)

1 bottle lasts approx. 50 patients.

**\$75** (50R205)



## Swipe Pads 96pk

Single use.

**\$25** (50R210)



## mirror magic® Organizer

Holds: 48 Swipe Pads and 1 - 6ml bottle of Anti-Fog Solution.

Store on operatory counter or in drawer.

Non-slip bottom for secure placement.

Note: box only, does not include mirror magic.

Dimensions: 3-3/4" x 3-5/8" x 1-7/8"

**\$16** (50R201)

S



**\$122 Value!**

**m**  
**w**  
mirror magic®  
anti-fog system



**Crystal HD**

## Clarity Collection

Includes:

Anti-Fog Solution 10pk

Swipe Pads 96pk

1 Crystal HD® Mouth Mirror

1 mirror magic® Organizer

**\$100** (50R220)

## The Perfect Pair

The superior visual acuity and 40% brighter Crystal HD® lens combined with the quick debris removal and Anti-Fog Solution of mirror magic® have proven to be the perfect pair. Increase chairside efficiency and speed up your process with this dynamic duo.



# Denture/Retainer Boxes



## Retainer Box

One piece construction.  
Solid locking mechanism.

### Retainer Box 1" Deep 12 pk

Dimensions: 3-1/8" x 2-3/4" x 1"

**\$13 (25R500) 12pk**

Color Code: J, N, O, P, R, S, Assorted

### Retainer Box 1-1/2" Deep 12 pk

Dimensions: 3-1/8" x 2-3/4" x 1-1/2"

**\$14 (25R550) 12pk**

Color Code: J, N, O, P, R, S, Assorted

S

High quality full color  
personalized imprinting  
now available!



## Denture Box

One piece construction.  
Solid locking mechanism.

### Denture Box 12 pk

Dimensions: 4-3/8" x 3-3/8" x 2"

**\$14.50 (30R800) 12pk**

Color Code: E, I, J, M, T, U, Assorted

S

## NEW! Personalized Imprinting

Our NEW and IMPROVED imprinting service is quick and easy. Just select a full color icon from our Fun or Professional sets or upload your own logo. Add three lines of text for your business name, phone number, website address, or your own unique message. It's as easy as that!

### Upload Your Own Logo Option

Actual size of imprinting area



### Denture Box Imprinting 24 pk

Dimensions: 4-3/8" x 3-3/8" x 2"

**\$41 (30P800) 24pk**

Color Code: A, E, I, J, M, T, U, Jewel Assorted

### Retainer Box 1" Deep 24 pk

Dimensions: 3-1/8" x 2-3/4" x 1"

**\$38 (25P500) 24pk**

Color Code: A, J, N, O, P, R, S, Vibrant Assorted

### Retainer Box 1-1/2" Deep 24 pk

Dimensions: 3-1/8" x 2-3/4" x 1-1/2"

**\$40 (25P550) 24pk**

Color Code: A, J, N, O, P, R, S, Vibrant Assorted

NOW  
AVAILABLE  
IN WHITE\*


\*Imprint Only



### Fun Imprinting Options



### Professional Imprinting Options



Visit us online to order!

## Saliva Ejector Accessories & MEGA Products



HVE Valve - Lever



Saliva Ejector Valve - Lever

**Hands free isolation!**  
**See page 22-23.**

### Saliva Ejector Valve (with Shut-Off)



Quick disconnect with swivel.  
Rubber tip with protective screen.  
Provides variable control from full on to full off.

### Saliva Ejector Valve with Lever On/Off Control

Dimensions: 3" x 5/8" x 3/4"

**\$35 (23E365) Gray Only**

### HVE Valve (Lever style)

Dimensions: 3" x 5/8" x 11/16"

**\$45 (23E300) Silver Only**

### Saliva Ejector Valve Lever Only

Dimensions: 5/8" x 3/16" x 5/8"

**\$10 (23E364) Gray Only**

### Saliva Ejector Replacement Tip

Dimensions: 5/8" x 5/8" x 1"

**\$8.50 (23E363) Gray Only**



### Saliva Ejector Screens



Single use.

### Saliva Ejector Screens

Dimensions: 3/8" Dia. x 3/8"

**\$10 (50Z952S) 25pk**

**\$70 (50Z952P) 250pk**



## MEGA Products

### MEGA Gracey

15-1/2" in length.  
Great teaching tool.  
A gift for any dental personnel.


**\$20 (50Z312) Mustard Only**



### MEGA Mirror

(Patient hand mirror).  
2-1/2" lens and 15-1/2" in length.  
Keeps children entertained.  
Great gift for any dental personnel.

**\$20 (50Z309) Color Code: N, P, Q, R, S**



**Fun gift ideas for any dental personnel. Keeps children entertained.**





## Water Wise & Steri-Soaker



### Water Wise Starter Kit

Holds 1 liter.  
Includes bottle, adaptor, lid, neon yellow and neon blue bubbles (water level indicator), and tube.

**\$33** (23Z100)

**S**



### Water Wise Replacement Bottle

Holds 1 liter.  
Includes bottle, blue lid, and neon yellow and neon blue bubbles (water level indicator).

**\$23** (23Z105)

**S**



### Water Wise vs. Traditional Water Bottle

	Water Wise	Water Bottle
Easy to clean	X	
Antimicrobial	X	
Easy-to-fill large mouth	X	
Quick release	X	
Level indicator color coded for: Sterile and distilled water, Medicament	X	
Strong	X	X
Easy to grip	X	
Sealable lid	X	
No threads to strip	X	
Pressure relief valve	X	

**We strive to make things easier around the dental office.**

## The Steri-Soaker

An efficient, high-volume germicide tray designed for maximum use of chemicals to clean/disinfect instruments. It is designed for easy cleaning and is large enough to accommodate ZIRC Steri-Containers, Cassettes, Bur Blocks, and Bur Guards.

### Steri-Soaker

Calendar on lid conveniently reminds staff when to change the solution; use color code rings to mark dates. Transfer basket's ribbed and perforated design allows quick drain-off and easy pick-up of instruments. Holds 3.78L of liquid; fill lines shown in quarts and liters.

Hinged fitted cover.

Safe-to-use chemicals: cleaners (enzymes, detergents), disinfectants (phenols)/cold sterilants for glutaraldehyde and chlorines.

Dimensions: 13-7/8" x 9-5/8" x 5-1/8"

**\$90** (50Z805) **Color Code: B, I**

**AM S**



**Designed for easy cleaning!**

# E-Z Access Shelf/Accessories

## E-Z Access Shelf

E-Z mounting to standard equipment post (1-3/4" - 2-1/4"); No tools needed.  
Will not damage post's finish.

E-Z adjustment to multiple positions.

Contoured edge to contain spills and materials.

Stain resistant finish; surface disinfect only.

Weight capacity of 25 lbs.

Workable tray area: 13-3/4" L x 10-1/4" W.


Standard E-Z Access Shelf as is accommodates 1-3/4" to 2-1/4" post.

Rotates 360°.

**\$125** (20Z420) Color Code: A, G, I



Provides easy  
access work space  
for auxiliary items.



## E-Z Access Shelf Adapters

### E-Z Access Shelf Adapter for Smaller Posts

Adapts to posts 1-3/8" to 2-3/4".

Dimensions: 2-1/8" x 1-7/8" x 3/4"

**\$5** (20Z422) Color Code: A, G, I

### 3-1/2" Bolt Adapter Kit for Larger Posts

Accommodates post sizes from 2-1/4" to 2-3/8".

Dimensions: 4" x 2" x 2"

**\$10** (20Z437) Color Code: A, G, I

### 4-1/2" Bolt Adapter Kit for Larger Posts

Accommodates post sizes from 2-3/8" to 2-3/4".

Dimensions: 4-1/2" x 2" x 2"

**\$10** (20Z438) Color Code: A, G, I



## Cup Holder-Snap On

Fits all standard cups.

No tools needed.

Dimensions: 3-1/8" x 2-1/2" x 1"

**\$11** (20Z436I)



## E-Z Access Shelf Surface Mount

Mount on horizontal flat surface only.

Dimensions: 3-3/4" x 3-3/4" x 2-7/8"

**\$17** (20Z424) Color Code: A, G, I



Our products are  
able to fit every  
need!



## 8" Extension Arm

Enables shelf to extend over patient.

Shelf has strength to hold 25 lbs. when connected to arm.

May be used with Wall Mount, Surface Mount, or directly to post.

Dimensions: 10" x 4" x 4-1/2"

**\$31** (20Z428) Color Code: A, G, I



## Wall or Cabinet Mount

Allows you to mount the E-Z Access Shelf onto a wall or the side of a cabinet. 90° angle of mount permits 180° movement of shelf on cabinet or next to wall.

Dimensions: 6-1/2" x 5-3/4" x 4-7/8"

**\$31** (20Z426) Color Code: A, G, I



# Cabinet Trays & Lab Pans

AM S

## Cabinet Trays



### No 16A – Hand Instruments (Regular Well)

Dimensions: 7-7/8" x 3-3/4" x 15/16"

\$8.50 (20Z201) Color Code: A, B



### No 16 – Hand Instruments (Deep Well)

Dimensions: 7-3/4" x 3-3/4" x 1"

\$8.50 (20Z202) Color Code: A, B



### No 17 – Clamps, Retainers, etc.

Dimensions: 7-7/8" x 3-3/4" x 3/4"

\$8.50 (20Z203) Color Code: A, B



### No 18 – Finishing Strips

Dimensions: 7-7/8" x 3-3/4" x 3/4"

\$8.50 (20Z204) Color Code: A, B



### No 19 – Pliers, Scissors, etc.

Dimensions: 7-7/8" x 3-3/4" x 3/4"

\$8.50 (20Z205) Color Code: A, B



### No 20 – Discs, Wheels, Stones

Dimensions: 7-7/8" x 3-3/4" x 3/4"

\$8.50 (20Z206) Color Code: A, B



Wall-Hanging

AM S

## Lab Pans

Organizes and transports impressions, models, dentures, and crown and bridge appliances. Stainless-steel spring clip holder. Color-coded.

### Lab Pan

Dimensions: 7-1/2" x 5-3/4" x 7"

\$6.50 (17Z101) Color Code: A, B, D, G, M

### Wall-Hanging Lab Pan

Same as standard version but with strip on top back for hanging on any vertical surface.

Saves on installing expensive shelving.

Dimensions: 7-3/4" x 5-3/4" x 7"

\$6.50 (17Z102) Color Code: A, B, D, G, M



## Wall Strip Lab Pan

Drilled and countersunk with screws and wall anchors. 18" Wall Strip.

Holds 3 Lab Pans.

Dimensions: 18" x 7/8" x 3/8"

\$9.50 (17Z104)

S

Save on installing expensive shelving.



## Dispensers



### Deluxe Cotton Roll Dispenser

Heavy-duty no-skid base.  
Dispenses up to three rolls at a time.  
Extra sturdy construction.  
Dimensions: 4-3/16" x 2-1/16" x 4"

**\$32 (28R800) Color Code: A, B**

S



### Cotton Roll Holder

Stainless pins embedded in weighted bottom.  
Neatly holds cotton rolls in an upright position.  
Hinged lid helps keep cotton rolls clean, prevents spilling.  
Dimensions: 3-5/8" x 3-1/4" x 2-3/8"

**\$30 (28R900) Color Code: A, B**

S



### Pellet Dispenser

Spring loaded – pushes pellets to the surface.

#### Dispenser for Small Pellets (sizes 3, 4, 5)

Dimensions: 1-7/8" x 1-7/8" x 2-1/16"

**\$31 (28R825) White Only**

#### Dispenser for Medium Pellets (sizes 0, 1, 2)

Dimensions: 1-7/8" x 1-7/8" x 2-1/16"

**\$31 (28R850) White Only**

S

### Sponge Dispensers

Spring action releases pads one at a time.  
Bottoms are weighted to prevent tipping.

#### 2" x 2" Sponge Dispenser

Dimensions: 2-1/2" x 2-1/2" x 2-7/8"


**\$31 (28R950) Color Code: A, B**

#### 4" x 4" Sponge Dispenser

Dimensions: 4-1/2" x 4-1/2" x 2-7/8"

**\$35 (28R975) Color Code: A, B**

S



### Matrix Dispenser

One roll of 3/8" x 12 yards Matrix clear tape 0.002 gauge.  
Allows Matrix to be contoured (u-shaped) as it dispenses.  
Stainless steel cutter allows for desired length.  
Disposable.

Dimensions: 2-1/8" x 2-1/4" x 5/8"

**\$26 (55Z200) Color Code: A, B**

S

# E-Z Storage/Handi-Hopper

## E-Z Storage Tub System

Hinges allow cover to stay open on its own!

Each Tub comes with one vertical divider and one half-depth adjustment divider.

Hinged cover comes in white or clear.

### E-Z Storage Organizer (with white cover)

Dimensions: 5-7/8" x 4-7/8" x 5-3/4"

Rack sold separately.

**\$34** (29R605)

### E-Z Storage Organizer (with clear cover)

Dimensions: 5-7/8" x 4-7/8" x 5-3/4"

Rack sold separately.

**\$34** (29R610)

### Replacement Hinges (pair)

Dimensions: 7/8" x 1/2" x 1/8"

**\$5** (29R627)

AM S

Hinges allow cover to stay open

Rack holding angled Tub



## E-Z Storage Single Rack

Holds storage organizer Tub.

Tub can be placed flat or at an angle for easy access of supplies.

Wall or cabinet mounted with optional mounting bracket.

Sturdy metal construction; plastic-coated.

White, non-slip feet.

Dimensions: 7" x 5" x 4-3/4"

**\$21** (29R650)

S

## E-Z Storage Mounting Bracket

White metal bracket for the E-Z Storage Tub single rack.

Pre-drilled holes for cabinet or wall mount.

Dimensions: 5-1/8" x 1" x 1/2"

**\$5** (29R651)

S



AM S

Mount this mini waste receptacle within arm's-reach of chairside for convenient disposal of contaminated waste.

## Handi-Hopper (with Mounting Bracket)

Small plastic waste receptacle for disposal of contaminated waste at chairside, keeping tray clean.

Slotted mounting bracket (included) for easy removal and replacement.

To use for infectious waste, order red (Biohazard) liners.

To use for other waste, order clear liners.

### Handi-Hopper

Dimensions: 7" x 4" 6-3/8"

**\$18.50** (20Z404) Color Code: G, I

### Handi-Hopper Bracket

Dimensions: 2-1/4" x 1/4" x 1"

**\$3** (20Z408) Color Code: G, I



## Handi-Hopper Liners (Single-Use)

Handi-Hopper Liners 100/Box (Red Biohazard)

**\$22** (20Z405) 100pk

Handi-Hopper Liners 100/Box (Clear)

**\$22** (20Z410) 100pk

AM S

## Waste-Drop Ring

Waste-Drop Ring for cabinet top.

Fits 5" hole.

Made of self-lubricating, easy-to-clean resin.

Dimensions: 5-5/8" Dia. x 3"

**\$14** (22Z100) Color Code: A, G, I



# Cuspidor/Dispos-A-Screen Traps & Needle Capper



**Reduce possibilities of infection and contamination.**

## Cuspidor Screens

Single use.

### Cuspidor Screens Large Universal – 2-1/2" (gray)

Adec, Ampco, Belmont, Dab Atlante, Dansereau, DCI Den-Tal-Ez, Forest, Knight, Marus, MDT, Proma, Ritter, Schein.

**\$18** (51Z205S) **Gray 25pk**

**\$71** (51Z205V) **Gray 144pk**

### Cuspidor Screens Small – 1-5/8" (white)

Adec, Cascade and Performer cuspidor bowls.

**\$18** (51Z200S) **White 25pk**

**\$71** (51Z200V) **White 144pk**

## Dispos-A-Screen Traps

Avoid contact with potentially hazardous material.  
Improves performance of suction system.  
Color coded by size. Single use.

**\$18** (25pk)

**\$70** (144pk)

### Blue – 2-1/8"

Adec models with black plastic canisters, plus newer models by Ampco, Belmont, DCI, Forest, Marus MDT, Schein.

**25pk** (51Z045S), **144pk** (51Z045V)

### Green – 1-7/8"

Ampco, Belmont, Biotec, Chayes, Dabi-Alante, Dansereau, DCI, Den-Tal-EZ, Dentech, Forest, Knight, Pelton & Crane, Proma and others.

**25pk** (51Z123S), **144pk** (51Z123V)

### Gray – 2-3/4"

Midwest 210 Models.

**25pk** (51Z070S), **144pk** (51Z070V)

*Left to right: Blue, Green, Gray, Yellow, Beige, Mauve and Coral*



### Yellow – 2-1/4"

Adec, model with metal canisters.

**25pk** (51Z040S), **144pk** (51Z040V)

### Beige – 2-7/8"

Den-Tal-EZ models with metal canister.

**25pk** (51Z060S), **144pk** (51Z060V)

### Mauve – 1-1/4"

Adec, Cascade, and Performer.

**25pk** (51Z080S), **144pk** (51Z080V)

### Coral – 2-5/8"

Pelton and Crane Spirit Unit with Plastic Canister.

**25pk** (51Z050S), **144pk** (51Z050V)



## Needle Capper

One-handed capping and recapping.

Syringe is for use; at convenient position in Needle Capper.

Base of device is secured to tabletop or tray with adhesive strip (included).

Needle is uncovered at same time syringe is withdrawn from Capper.

After use, syringe is replaced for recapping.

Needle Capper available in Standard and Long sizes.

Either size tube will fit standard base.

### Standard Needle Capper

(base, one tube, adhesive tape)

Dimensions: 2-3/16" x 1-1/4" x 1-1/4"

**\$26** (27R700) **White Only**

### Long Needle Capper

(base, one tube, adhesive tape)

Dimensions: 2-1/2" x 1-1/4" x 1-3/4"

**\$26** (27R704) **White Only**





# Mighty Mixer Bowls, Spatulas & Organizers



Single-use  
disposable mixing  
bowls.



## Mighty Mixer Bowls 36pk (White Only)

Disposable, eliminating cleanup time.  
Fits securely into mixing bowl without slippage.  
Saves wear on mixing bowl.  
Free-standing for use with or without mixing bowl.  
Capable for use with electric Alginate mixers.  
Single-Use.  
Dimensions: 4-7/8" Dia. x 3-1/8"

**\$21 (50Z505) 36pk**



## Mighty Mixer Bowls 36pk (Assorted)

Same features as the white Mighty Mixer Bowls.  
Dimensions: 4-7/8" Dia. x 3-1/8"  
9 of each color: J, N, P, and R  
Single-Use.

**\$21 (50Z505-ASTD) 36pk**

## Mighty Mixer Spatula 4pk (Assorted)


Same features as the white Mighty Mixer Spatula.  
Dimensions: 7-5/16" x 3/16" x 1-1/4"  
1 of each color: J, N, P, and R

**\$20 (50Z501-ASTD) 4pk**



## Mixing Sticks 100pk

Single use.  
Dimensions: 4-1/2" long  
**\$21 (50Z521) 100pk**



## Cement Spatula

White plastic, double-end spatula.  
For fillings, cement, and pastes.  
Self-lubricating plastic for easy cleaning.  
Dimensions: 6-3/8" x 1/4" x 3/8"

**\$6 ea. (50Z502)**



## Delrin Spatula 10pk

Flexible flat blade.  
Non-stick.  
Resistant to chemicals.  
Great for mixing cements, acrylics, etc.  
Dimensions: 4-5/8" x 3/16" Dia.

**\$10 (50Z523) 10pk**



## Mighty Mixer Spatula

Self-lubricating for easy cleaning.  
White plastic spatula.  
Broad, medium flex blade.  
Ideal for alginates and plaster.  
Dimensions: 7-5/16" x 3/16" x 1-1/4"

**\$6 ea. (50Z501) White Only**



No more  
confusion!



## Implant Organizer

Quick and easy recognition. Individual compartments for implants.  
Shape of the organizer matches the arch of the mouth for implant placement.  
Separate compartments prevent implant appliances from mixing together.  
Lubricated plastic enables tooth number to be written in the designated area.  
Cover helps keep materials secure.  
Dimensions: 6-7/8" x 4-5/8" x 5/8"

**\$25 (20Z490B) Blue Only**



## Restorative Organizer

Compartments shaped like the arch of the mouth for veneers, inlays, composite, and onlays.  
Disposable Mixing Wells included for bonding agent and cement.  
Accessory area for brushes and materials.  
Lubricated plastic enables tooth number to be written in the designated area.  
Light protective cover can be used as a light shield and to prevent materials from setting early.  
Includes 25pk of Mixing Wells.  
Dimensions: 6-7/8" x 4-5/8" x 5/8"

**\$30 (20Z495R) Neon Purple Only**



## Mixing Well (Single-Use) 144pk

Dimensions: 2-7/8" x 1-1/8" x 3/16"

**\$20 (20Z496V) 144pk**



\*All US pricing subject to change

Contact us today to get started!  
800.328.3899 | [sales@zirc.com](mailto:sales@zirc.com)