

ISOLATION360™

POWER PACK



Serves as a barrier to swallowing or aspirating objects, allowing dentists to practice safely and confidently

Ideal for Crown & Bridge, Endodontics, Extractions, Implants, Inlays/Onlays, Orthodontics, Veneers

- Patient protection
- Easy placement & removal
- 3 convenient sizes (Small, Medium & Large)
- Single Use

SE CUSHION™

Ideal for Prophyls, Endodontics in conjunction with Insti-Dam®, Whitening, Orthodontics

- Prevents tissue grabbing and tissue irritation
- Increases patient comfort
- Improves range of suction
- Quiets the noise of suction
- Single Use

INSTI-DAM®

Ideal for Endodontics

- Built-in flexible frame allows easy access for x-rays without removal
- Pre-punched hole makes positioning easy and helps eliminate tearing (additional holes may be punched)
- Compact design fits outside patient's lips
- Easy-to-use, and comfortable for patient
- Latex Free
- Single use and pre-assembled



Ideal for Scaling & Root Planing & Ultrasonic Scaling

- Holds suction tip in place
- Protects and retracts cheek
- Patient can close around suction with ease
- Easily slides on the saliva ejector
- Soft material for patient comfort
- Single Use

LINGUA-FIX™

Ideal for Scaling & Root Planing, Digital Imaging, Ultrasonic Scaling, Sealants

- Provides a dry field
- Retracts and protects tongue
- Protects tongue from injuries
- Easy to place & remove
- Comfortable, no sharp edges or corners
- Large suction surface avoid blockage
- Single Use and pre-assembled



Ideal for Prep of Crowns, Bridges, Inlays, Onlays and Veneers, Composite / Amalgam, Sealants


- Gain profitable time in your day
- Isolate, retract and evacuate during almost any procedure Hands-Free!
- Efficient suction creates a dry and ideal working environment
- Trimmable or a custom fit
- Connects directly to HVE valve
- Built-in bite block
- Single Use and pre-assembled

INCLUDES

- 5 Lingua-Fix™
- 10 Pink Petal™
- 20 SE Cushion™
- 9 Airway Armor™
(3 Small, 3 Medium & 3 Large)
- 10 Mr.Thirsty One-Step®
(5 Pedo/Small & 5 Med/Large)
- 10 Insti-Dam®
(5 Original & 5 Relaxed Fit)

\$100 | Isolation360™ Power Pack
(50Z950)

50Z950LT-1_REVA_0918



ISOLATION360™
your total isolation partner

7 AWESOME ISOLATION DEVICES, ONE HAPPY YOU

This is your complete Power Pack for all your isolation needs!

The ultimate goal is patient safety, comfort and practice efficiency. The pack includes 7 of our innovative products that are made in the U.S.A.

Zirc's Isolation 360™ has you covered!



Zirc
Because Time Is Everything

www.zirc.com | 800.328.3899