



PRODUCT GUIDE 2019



## WHY ZIRC

2

## COLOR METHOD™

4 - 19

Why the Color Method™?	4
Color Chart	4
Zobe	5
Tub Kits	7
Tub Accessories	8
Trays   Covers	10
Multi-Mod Racks	12
Cassettes	13
Steri-Containers   Mats	15
Bur & Endo Guards/Blocks	16
E-Z ID Rings/Tape	19

## ISOLATION

20 - 27

Isolation360™ Power Pack	20
Airway Armor®	22
Pink Petal®	23
mr.thirsty® one-step	24
Insti-Dam®   SE Cushion   Lingua-Fix®	26

## MOUTH MIRRORS

28 - 33

Crystal HD® Mirrors	28
mirror magic® Anti-Fog Solution	32

## OPERATORY/MISC

34 - 44

CAD/CAM Block Locker™	34
-----------------------	----

# MAKING DENTAL TEAMS HAPPIER AND MORE EFFICIENT




## Zirc History - 1967

James Campion, our CEO, founded Zirc in 1967 and his daughter Nicolle, who joined the team in 2000, is the company's president. In 2017 we celebrated 50 years of offering our customers the highest quality products at the best value. Dental technology and processes have seen many changes over the last 50 years, but at Zirc, our commitment to providing time-saving, innovative dental products continues to be our top priority. We look forward to serving our customers for years to come!



## Color Method™

Zirc has a brilliant solution for organizing the dental practice by the method of color-coding. We have rebranded our color-coding system to the Color Method™, to better educate the concept of color-coding in the practice. The method gives step-by-step instructions to achieve an efficient, productive, and compliant work place...all while remaining simple and stress-free.



## Made in the USA

Our headquarters, located in Buffalo, Minnesota, houses our engineering, tooling, manufacturing, distributing, and customer service operations. We're proud to say that our state-of-the-art, high-quality products are made in the U.S.A.



## Antimicrobial


Antimicrobial compounds are intended to protect parts against staining, odors, and degradation due to fungi, mold, yeasts, mildew, algae, and other microbes. Zirc incorporates antimicrobial compounds during the manufacturing process. Unlike a coating, this prevents fading or wearing away, giving you long-lasting protection. The incorporation of antimicrobial compounds in our plastic products is not intended to protect users or others against food-borne or disease-causing bacteria. For more information, please reference EPA Registration Notice 2000-1.



## Zirc Perks

Members receive FREE samples, money-saving promotions, tips to get and stay organized, product evaluation invites, and more! Sign up to become a Zirc Perks Member via [www.zirc.com](http://www.zirc.com) or call us at 800-682-6636 (available in U.S.A. only).





## ASPIRATION & INGESTION PREVENTION

### Airway Armor®

Airway Armor® prevents foreign objects and debris from going down the throat during procedures, giving yourself a peace of mind and your patients the protection they deserve. It is easy to place, tear resistant, and comfortable to the patient, all while giving you full mouth accessibility. See page 22.

Find the products  
that fit your needs!



## CONTROL BLOCK CHAOS

### CAD/CAM Block Locker™

Our CAD/CAM Block Locker™ is a stylish, modern cabinet designed to help organize the chaos of milling blocks. Utilize the online template to label each internal organizer tray by shade, color and opacity...creating a seamless way to check inventory, stock and keep blocks organized with ease. See page 34.



### Safe-Lok® Cover Tall

Maintain compliance with OSHA and CDC with our locking tray cover! Now offered in a tall version, you will be able to fit our Performance Tool Organizer® and almost any stainless steel cassette! See page 11.

## What's New?

Continuous,  
hands-free  
suction!



## GIVE A HAND TO HYGIENE

### Pink Petal®

Provides continuous, comfortable, hands-free suction with no additional attachments required. By placing buccal side, the Pink Petal® confidently holds the saliva ejector right where you need it. See page 23.

## STRUGGLING WITH ISOLATION?

### Isolation360™ Power Pack

Zirc's Isolation360™ has you covered! The pack contains 7 different isolation devices, allowing you to become familiar with each one. Also included inside the pack, is an Ideal usage chart with product recommendations for different procedures.

Once you use all the products, order more of what works best for your practice! See page 20.



For convenience, safety and efficiency!

## COMING SOON!

### Pouch-it™

Zirc will be releasing a line of custom-fit pouches that will pair with our instrument holders, eliminating the need of guessing what size to get! See page 13.





**Organization is the goal.  
Color is the method.  
Zirc is your solution.**

## ORGANIZE WITH EASE Why the Color Method™?

**COLOR METHOD™**  
THE KEY TO ORGANIZATION



Zirc has the tools you need to manage your materials and easily store, access, and clean your instruments. Staff, and patients alike, will appreciate a streamlined, stress-free day, from set-up to sterilization. Our exclusive Color Method™ will help you get organized, increase efficiency and eliminate leaving chairside to search for materials.

**Color coding has been proven to save time for each procedure, allowing you to see more patients in a day.**

### Simplified Practice

Once you implement the Color Method™, your workflow will become a systematic routine between the sterilization and treatment room, creating maximum efficiency!

### Benefits

- Set-ups are organized, eliminating the need to find missing materials or instruments during a procedure
- Makes re-stocking dental materials more efficient by decreasing the amount of inventory, and reducing the amount of expired materials
- Simplified training of new staff
- Reduces the risk of an exposure incident by containing dirty instruments in covered a container, satisfying the CDC and OSHA standards
- Visually pleasing to patients versus stainless steel

*For more information, visit us online to download our Color Method™ brochure.*

### Color Chart

With 16 product colors, there is guarenteed to be a combination thats matches your practice!



# Want help implementing the Color Method™?



**Z**IRC'S **O**RGANIZATIONAL **B**OX OF **E**FFICIENCY is a complementary resource that assists you in implementing the Color Method™. It contains over 30 of our top-selling organizational products that will help you determine what's needed for specific procedures within your practice.

After a short consultation with a Zirc Efficiency Expert, we will send you the ZIBE. You'll be able to have a hands-on experience with the products, while figuring out what products work best for your practice. Once you are ready to proceed, we will help customize your order! Then use the prepaid return label to send us the ZIBE back.

*Note: The products contained in the ZIBE are for exploration of products only, they are not intended for clinical/procedure use. Only available in the U.S.A.*

## ORDER TODAY!

visit [www.zirc.com/zobe](http://www.zirc.com/zobe)  
or call 800.328.3899



what's in  
the box?!







Syringe Composite Tub



Capsule Composite Tub



Endo Tub



Crown Prep Tub



Hygiene Tub



Implant Tub

## SMALL CHANGES Big Results

**Cluttered drawers = stress and chaos. Stay efficient with color-coded organization.**

Start by removing everything from your operatory drawers. Place in a large box and set aside. Restock your operatory drawers with staple supplies such as: bibs, saliva ejectors, cotton products, etc. (basically, anything you use in EVERY procedure).

Here's where the Tubs come into play. You need a minimum of one Tub for each procedure. Each Tub should be a different color so that you can quickly identify procedures by color. Stock each Tub with all materials you will need for that procedure.

Store ALL Tubs in a central area (such as sterilization). Take the Tub into the treatment room when a procedure is performed. Pull out what you need at chairside. Disinfect and re-store.

**Reduce your stress and have fun  
with color organization.**





## MAKE IT SIMPLE WITH AN All-In-One Kit

Our kits are designed to store and organize materials and keep your inventory under control. In addition, they are easily portable for transportation between procedures and your centralized location for sterilization, disinfecting, and restocking.


### Complete Tub

Includes: Procedure Tub, Tub Liner, Tub Dividers, Cup/Cover Single, Cup/Cover Double, Divided Slide Tray, Safe-Lok Cover for Tub®.

Dimensions: 12-5/8" x 11-1/4" x 4"

**\$68 (20Z455) Available in all colors**

AM S



Use: Any Procedure

## Composite Kits


Store and organize composite materials with our convenient Composite Kits. When stored in a Tub, contents become portable and can be transported to a central sterilization area for restocking. Organizational accessories help you keep better inventory control with Clip-On Holders that allow you to separate composite for use during procedures. The Long Tub Cup and Slide Tray provide additional organizational space. The Shade Label Sheet provides self-stick labels preprinted with composite shade.

### Syringe Composite Kit

Syringe holder allows you to store up to 40 syringe tubes. Kit includes: Syringe 20-Unit Stand, 1 Long Tub Cup/Cover, 1 Divided Slide Tray, 1 Syringe Clip-On Holder, Shade Label Sheet 2pk.

**\$52 (20Z489)**

AM S



### Capsule Composite Kit

Long Tub Cups provide a total of 20 compartments for capsule storage. Includes: Tub Liner, 4 Long Tub Cup/Cover, 1 Divided Slide Tray, 1 Capsule Clip-On Holder, Shade Label Sheet 2pk.

**\$52 (20Z479)**

AM S





## Build Your Own Tub



### Procedure Tub

Allows you to “build” your own setup.  
Dimensions: 12-1/4" x 10-7/8" x 2-3/4"

**\$19.50** (20Z463) Available in all colors



### Safe-Lok® Tub Cover

Locking cover for securing materials.  
Keep materials protected during storage and transportation.  
Dimensions: 12-1/4" x 10-7/8" x 2-3/4"

**\$19** (20Z459)



### Tub Liner

Foam liner that sits in bottom of Tub.  
Tub not included. Replace as needed.  
Dimensions: 9-5/8" x 8-5/8" x 1/8"

**\$6** (20Z462)



### Tub Box

Designed to hold/organize “hard-to-find” supplies.  
Non-slip bottom for secure placement.  
Dimensions: 3-3/4" x 3-5/8" x 1-7/8"


**\$15** (20Z474) Gray Only



### Tub Insert

10 compartments for easy organization.  
Sits on rim of Tub, allowing for storage underneath.  
Dimensions: 10-3/4" x 9-11/16" x 1"

**\$15** (20Z465) Gray Only



### Tub Dividers

Creates up to 24 compartments.  
May be cut or trimmed to create larger compartments.  
Dimensions: 10-1/16" x 9-1/16" x 1-5/16"

**\$12.50** (20Z461) Gray Only



Use color-coded tape on Tub Cups and Tub Box to identify procedural contents. (Learn more on page 19)



### Tub Cups with Covers

Fill Tub Cups with brushes, clamps, cotton pellets, prophylaxis angle tips, prophylaxis paste, pumice, syringe tips, implant tools, etc. Can be used in procedure tubs, drawers, on trays, etc. for small, easy-to-lose items. Covers are transparent.



### Single Tub Cup with Cover

Dimensions: 2-1/4" x 1-5/8" x 1-1/8"

**\$7.50** (20Z471) Gray Only

### Double Tub Cup with Cover

Dimensions: 3-3/8" x 2-1/4" x 1-1/8"

**\$7.50** (20Z472) Gray Only

### Long Tub Cup with Cover

Dimensions: 8-3/8" x 2-5/16" x 1-1/8"

**\$13.50** (20Z473) Gray Only

# Tub Accessories



## Divided Slide Tray


Dimensions: 9-3/4" x 5-3/8" x 1-1/16"

**\$12** (20Z466) Clear

## Flat Slide Tray

Dimensions: 9-3/4" x 5-3/8" x 1-1/16"

**\$12** (20Z467) Clear



## Shade Label Sheets 2pk

Preprinted labels for marking composite materials within the Tub organizers.

Dimensions: 7" x 4-3/8" x 1/64"


**\$5** (20Z477)

## Slide Trays

Sits on rim of Tub and allows for storage underneath.



*Check out our Implant and Restorative Organizers on page 44.*



## Syringe 20-Unit Stand

Stands up for easy access to composite.

Lays flat to fit inside Tub or drawer.

Holds up to 20 syringes.

Has space for 2 Syringe 10-Unit Add-Ons (Stand can hold up to 40 syringes).

Dimensions: 8-3/4" x 4-3/4" x 1-9/16"

**\$25** (20Z487)



## Syringe 10-Unit Add-on

An Add-on for the Syringe 20-Unit Stand.

Gives additional space for syringe storage.

Holds 10 syringes.

Dimensions: 8-9/16" x 2-5/16" x 3/4"

**\$10** (20Z488)



Syringe 20-Unit Stand

*The Syringe 10-Unit Add-on may be added for storage of up to 40 syringe tubes.*



## Syringe Clip-On Holder

Allows quick access to composite during procedure.

Conveniently clips on any side of the Tub or B-Lok Tray.

Holds 15 syringes and 3 applicator tips.

Dimensions: 8" x 2-7/16" x 3/4"

**\$10** (20Z486)



## Capsule Clip-On Holder

Allows quick access to composite during procedure.

Conveniently clips on any side of the Tub or B-Lok Tray.

Holds 8 capsules, 4 syringes, and capsule/gun dispenser.

Dimensions: 8" x 2" x 9/16"

**\$10** (20Z476)





## Trays & Covers



### Trays

Our exclusive Safe-Lok® tray cover provides an excellent method of infection control and safety that allows for two benefits other trays do not:

- 1** Prevents any contaminated items from falling off of the used tray during transportation from operator to central sterilization. The contaminated items stay enclosed until you are ready for the cleaning and sterilization process.
- 2** Set up your procedure trays ahead of time with all instrumentation needed. The Safe-Lok® cover securely locks to the tray, not only keeping items clean, but also allowing for safe transportation to chairside. The cover is removed once the patient is in the chair and the procedure begins.

### Are you in compliance?

#### OSHA

Immediately or as soon as possible after use, contaminated reusable sharps shall be placed in puncture-resistant, leak-proof and labeled or color-coded containers until properly reprocessed. (To avoid spillage of contents, it is suggested that they be covered and secured prior to moving.)

#### CDC

Dirty Instruments must be "contained" prior to transporting them. Containment can be defined as **instrument cassettes** or cages in which dental care items are cleaned, sterilized, and stored until point of use. Sharp instruments should not be carried openly to the instrument-processing area. Percutaneous injuries may occur to other dental healthcare personnel (DHCP) or patients in hallways leading to the instrument-processing area.

**The Safe-Lok® Tray Cover helps prevent cross contamination and provide safe transportation back to the sterilization area.**

# Trays & Covers



## B-Lok Flat Tray

Wide edges for secure grip.  
Non-slip feet for stability.  
Self-lubricating, high-gloss surface helps prevent staining.  
Dimensions: 13-3/8" x 9-5/8" x 7/8"

**\$13 (20Z401) Available in all colors**



## B-Size Safe-Lok® Tray Cover

Fits B-Lok Flat and Divided Trays only.  
Side locking mechanism securely locks cover to tray.  
Dimensions: 13-7/8" x 9-7/8" x 3/4"

**\$18 (20Z445)**



**NEW!**



## B-Lok Divided Tray

Holds 12 hand instruments.  
Bur holder for 9 FG, RA, or HP burs.  
Comes with Insti-Grip.  
Non-slip feet for stability.  
Self-lubricating, high-gloss surface helps prevent staining.  
Holder for Syringe Tip or Saliva Ejector.  
Dimensions: 13-3/8" x 9-5/8" x 7/8"

**\$15 (20Z451) Available in all colors**



## B-Size Safe-Lok® Cover Tall


The Tall size accommodates our Performance Tool Organizer and almost any stainless steel cassettes.  
Fits B-Lok Flat and Divided Trays only.  
Side locking mechanism securely locks cover to tray.  
Dimensions: 13-7/8" x 9-7/8" x 1-3/16"

**\$18 (20Z446)**



*Our Performance Tool Organizer™ fits on Zirc Trays with the Safe-Lok® Tall Cover! See page 13 for E-Z Jett® Cassettes.*

## FOR B-LOK DIVIDED TRAY



## Insti-Grip

To be used with the B-Lok Divided Tray only.  
Holds instruments quietly and safely in place.  
Dimensions: 6-5/8" x 1/4" x 3/8"

**\$2 (20Z454)**

## B-Size Tray Cover (Non-Locking)

Dimensions: 14" x 10-1/8" x 5/16"

**\$18 (20Z441)**



## Mini Tray

Wide edges for secure grip.  
Non-slip feet for stability.  
Self-lubricating, high-gloss surface helps prevent staining.  
Mini Trays are great for hygiene/exam setup.  
Dimensions: 9-3/8" x 6-3/8" x 7/8"

**\$12 (20Z101) Available in all colors**



## Mini-Lok Tray Cover

Side locking mechanism securely locks cover to tray.  
Dimensions: 9-7/8" x 6-5/8" x 1"

**\$15.50 (20Z103)**



## Mini Tray Cover (Non-Locking)

Dimensions: 9-7/8" x 6-5/8" x 1"

**\$15 (20Z102)**





# Multi-Mod Racks

## SAVE ON SPACE & Organize Your Office

Tub, tray, and cassette racks are very effective additions to a dental practice. They save on counter space and several tubs or trays with covers can fit into one rack. These racks have non-slip feet so they are very stable on any counter or cabinet surface. With a chrome-plated finish, these racks can fit any dental practice design.

### Multi-Mod Rack 8-Place

Holds 8 Trays or 4 Tubs with covers.  
Adjustable to hold Trays, Tubs, and Cassettes.  
Sturdy metal construction; chrome-plated.  
Non-slip feet.

**\$94** (21Z410)

S

### Multi-Mod Rack 6-Place

Holds 6 Trays or 3 Tubs with covers.  
Adjustable to hold Trays, Tubs, and Cassettes.  
Sturdy metal construction; chrome-plated.  
Non-slip feet.


**\$94** (21Z105)

S

Multi-Mod Rack (8-Place)



Multi-Mod Rack (6-Place)



**Adjustable  
to fit your needs!**



*Can adjust to hold  
Mini Trays and Cassettes*

## WHY USE Zirc's Cassettes?

A low-cost investment with a quick return! Individual instrument handling takes time away from treating patients and increases overall processing time. Zirc's cassettes will prolong your instrument's life and save you time without breaking the bank.

In a 2014 Implementation of Cassettes Case Report, Marie T. Fluent, DDS, writes how the dental office benefited from the transition from no cassettes to using them with every procedure: "Drs. Dale and Thomas were both grateful that they chose to implement ... an Instrument Management System for the organization, processing, and handling of their instruments within their practice. They agreed that cassettes improved all aspects of instrument management... Initially, their opinions were varied; however, at the conclusion of the one-year anniversary, no clinical staff could imagine processing instruments in any other fashion. Cassette usage had become organized, efficient, safe, and routine." Read the article in its entirety at: [www.dentaleconomics.com](http://www.dentaleconomics.com).

- Designed for superior water flow and cleaning - open top, round holes, and radius corners
- Keeps instrument setups together through sterilization and procedure processes
- Protects instruments from extensive damage and prolongs their life
- Non-slip feet for stability and safety
- AUTOCLAVABLE
- Available in 16 colors for color coding by procedure

Keep your workplace  
organized and confusion  
to a minimum!



Learn more at  
[www.zirc.com/instrumentcassettes](http://www.zirc.com/instrumentcassettes)

## Sterilization Pouch Size

Cassette	Pouch Size
E-Z Jett® 5-Place	5-1/4" x 12"
E-Z Jett® 8-Place	7-1/8" x 13-13/16"
E-Z Jett® 8-Place Slim	7-1/8" x 13-13/16"
E-Z Jett® 10-Place	10" x 15"

coming soon!

### Pouch-It™

Custom-fit sterilization pouches that will eliminate the guess work of what size pouch to use with each Cassette

## Sizing Chart for E-Z Jett® Cassette Sterilization

Autoclave Manufacturer	Model	Chamber Dimensions L x W x D	E-Z Jett® 5-Place	E-Z Jett® 8-Place	E-Z Jett® 8-Place Slim	E-Z Jett® 10-Place
Midmark	M11	11" x 18"	12	8	8	4
	M9	9" x 15"	6	4	4	2
Pelton & Crane	Magnaclave	15-1/4" x 15" x 26"	26	24	24	12
	Delta XL 10	9-13/16" x 17-1/2"	12	9	9	4
	DeltaXL8	8-1/4" x 14"	8	4	4	2
SciCan	Bravo17	10" x 13-1/2"	9	6	6	3
	Bravo21	10" x 17-1/2"	12	9	9	6
Tuttanuer	3870	15" x 27"	20	12	12	10
	EZ10	10" x 18"	12	9	9	6
	2540M	10" x 18"	12	9	9	6
	EZ9	9" x 18"	8	6	6	4



## AUTOCLAVABLE Instrument Holders




Syringe Tip Holder

**Syringe Tip Holder 5pk**  
Dimensions: 1" x 3/8" x 3/8"  
**\$5 (50Z926)**

### E-Z Jett® Cassette 5-place

Holds 5 instruments.  
Great for a basic exam setup!  
Syringe Tip Holder included.  
Covers stay attached to the cassette and fold back easily.  
Fits Statim Sterilizer.  
Dimensions: 8" x 3-1/8" x 1-3/16"


**\$28 (50Z925) Available in all colors**



### Instrument Mat

Holds 10 hand instruments.  
Reverse side holds 4 surgical instruments.  
Non-slip, soft material.  
Small size allows easy access to instruments.  
Organizes instruments in order of procedure use.  
Can be cut to fit Mini Tray.  
Dimensions: 6-3/4" x 2" x 3/8"


**\$14.50 (20Z464) Available in all colors**



### E-Z Jett® Cassette 8-place

Holds 8 instruments.  
Covers stay attached to the cassette and fold back easily.  
Syringe Tip Holder included.  
Dimensions: 7-15/16" x 4-7/16" x 1-1/8"

**\$38 (50Z927) Available in all colors**



Instrument Strap

### Mini Mat


Holds 8 instruments.  
Organizes instruments at chairside and for sterilization.  
Non-slip feet for stability.  
Dimensions: 3-7/8" x 3-1/2" x 1/2"

**\$13 (20Z480) Available in all colors**

### Instrument Strap 2pk

Locking strap secures instruments.  
Dimensions: 3-5/8" x 1/2" x 5/16"

**\$6/pair (50Z917) Gray Only**




*Size comparison with the  
8-Place E-Z Jett® Cassette.*



### E-Z Jett® Slim 8-place

Holds 8 instruments.  
Covers stay attached to the cassette and fold back easily.  
Compact size. Fits the Statim.  
Stackable in Ultrasonic.  
Syringe Tip Holder included.  
Dimensions: 7-1/16" x 4-1/16" x 1-5/16"

**\$33 (50Z911) Available in all colors**



### E-Z Jett® Cassette 10-place

Holds 10 instruments.  
Accessory compartment in the middle for extra storage.  
Covers stay attached to the cassette and fold back easily.  
Dimensions: 7-15/16" x 7" x 1-1/8"

**\$43 (50Z930) Available in all colors**



# Steri-Containers



Closed

## Performance Tool Organizer

A deeper, roomier cassette made for large and specialty instruments. Design your own layout with our customizable, removable feet (also sold separately). Designed for superior water flow for superb cleaning. Compact styling to fit in most sterilizers. Covers stay attached to the cassette and fold back easily. Ideal for:

- Surgical elevators and forceps
- Orthodontic pliers
- Hemostats
- Scissors
- Any instrument that does not fit in your regular cassette



Open

### Performance Tool Organizer

Covers stay attached to the cassette and fold back easily.  
Dimensions: 8-3/16" x 4-1/4" x 1-7/32"

**\$56** (50Z940) Available in all colors



### Performance Tool Organizer Replacement Feet

Two strips of 10.

**\$7** (50Z944) Gray Only



**Holds your large and specialty tools**

*E-Z Jett® Cassettes fit on Zirc Trays with locking cover and the Performance Tool Organizer® fits our new Safe-Lok® Cover Tall! Check out our Trays & Safe-Lok® Covers on page 10.*



## Steri-Container

Holds 10-14 instruments.  
Keep instrument setups together through sterilization processes.  
Dimensions: 8-1/8" x 1-7/8" x 1-7/8"

**\$11.50** (50Z900) Available in all colors



## Compact Steri-Container

Holds 8-12 instruments.  
Easy open latch, just push down.  
Keep instrument setups together through sterilization processes.  
Statim 2000 holds 7 Compact Steri-Containers.  
Dimensions: 7-1/8" x 1-1/2" x 1-1/2"





# Bur Guards



## PROTECT YOUR BURS Steri-Bur Guards




### Steri-Bur Guards and Adapters

This gadget has it all! Not only does it protect your burs during cleaning and sterilization, but it can also be positioned at three different angles to give you easy access at chairside.

#### Steri-Bur Guard

- Fits FG, RA, and HP Burs
- Adjustable height for burs of different lengths (12-Hole Tall Bur Guard not adjustable)
- Holds burs in position for maximum cleaning effectiveness
- Drain holes help prevent rusting of burs



#### 12-Hole Bur Guard

Accommodates burs from 12mm-30mm.  
Dimensions: 2-7/8" x 3/8" x 1-3/8"

**\$15.50 (50Z406) Available in all colors**



#### 22-Hole Bur Guard

Accommodates burs from 12mm-30mm.  
Dimensions: 4-7/8" x 7/16" x 1-7/16"

**\$16 (50Z410) Available in all colors**



#### 12-Hole Tall Bur Guard

Accommodates burs from 30.9mm-37.1mm.  
Dimensions: 2-7/8" x 3/8" x 1-5/8"  
Ideal for nose-cone lab burs



**\$16 (50Z408) Gray Only (I)**

*A Short Bur Adapter is included to stop any shorter burs from falling out.*

#### Universal Bur Adapter 5pk

Can be used with 12- or 22-Hole Bur Guard.  
Enables long and short burs to be used in the same bur guard without short burs falling out.

**\$5 (50Z412)**



# Endo Guards & Holders




## Endo Guards

A unique delivery system to organize, store, and sterilize endodontic files and reamers. The Endo Guard swivels 360° with stops at every 45° for operator's choice. Large slotting provides complete visibility and draining during cleaning.

### Steri-Endo Guard

- Holds most lengths for staging, cleaning, and sterilizing files and reamers
- Millimeter gauge measures and positions rubber stopper
- Counter on locking bar tracks the number of sterilizations
- Right- or left-hand use; optional angling
- Available in Hand or Engine styles
- Single-use foam insert cleans/holds files (packaged with 12)



### Steri-Endo Guard Engine

Organizes up to 11 engine and 5 hand files/reamers.  
Dimensions: 5-3/8" x 3/8" x 2-3/16"

**\$31 (50Z452) Available in all colors**

### Steri-Endo Guard Hand

Organizes up to 16 hand files/reamers.  
Dimensions: 5-3/8" x 3/8" x 2-3/16"

**\$31 (50Z450) Available in all colors**

### Foam Inserts

Single use.

**\$31 (50Z451) 48pk**

**\$62 (50Z451V) 144pk**



*Unique organizer rests on your finger.*

### Endo Assist

Built-in millimeter gauge measures and positions rubber stopper. For right- or left-hand use. Small and large rings are included to accommodate most finger sizes. Flat bottom enables the Endo Assist to stand freely on its own. Single-use foam insert cleans/holds files (packaged with 12). Dimensions: 1-1/2" x 1-3/4" x 1-5/8"


**\$28 (50Z460) Available in all colors**



### Assist Stand

Dense foam cleans bio-debris faster and more easily. Holds Endo Instruments in an upright position. Square shape fits great in procedure tub. Thicker foam prevents solution from leaking. Rounded corners for easy cleaning. Dimensions: 1-7/8" x 1-7/8" x 1-3/8"

**\$16.50 (50Z462) Black Only**



### Assist Stand Foam

Single use.

Dimensions: 10-1/2" x 7-3/4" x 7/8"

**\$31 (50Z466) 48pk**

Dimensions: 10-1/2" x 7-3/4" x 7/8"

**\$62 (50Z466V) 144pk**



## Bur Blocks



### Magnetic Bur Blocks

- Fits FG, RA, and HP Burs
- Magnetic base holds burs securely
- Non-slip feet for stability
- Drain holes help prevent rusting of burs



#### 8-Hole Bur Block

Covers sold separately.  
Dimensions: 1-1/2" x 1-1/2" x 1/2"

**\$14.50** (50Z400) Available in all colors



#### 14-Hole Bur Block

Covers sold separately.  
Dimensions: 2-15/16" x 1-1/2" x 1/2"

**\$15.50** (50Z403) Available in all colors

Recessed area to hold used burs!



### Bur Block Covers

- Covers available to maintain cleanliness of burs during storage
- 1" covers for FG and RA burs; 2" cover (14-hole only) for lab and surgical burs
- Covers are not autoclavable
- Surface disinfect only



#### 8-Hole Bur Block Covers 1" High - Standard

Dimensions: 1-5/8" x 1-5/8" x 1-1/16"  
**\$7.50** (50Z401)



#### 14-Hole Bur Block Covers 1" High - Standard

Dimensions: 3-1/8" x 1-5/8" x 1-1/16"  
**\$7.50** (50Z404)



#### 14-Hole Bur Block Covers 2" High - HP Lab Burs

Dimensions: 3-1/8" x 1-5/8" x 2-1/4"  
**\$7.50** (50Z405)

# Rings & Tape

## How to effectively identify your instrument setup.

**1 color for procedure:**  
Vibrant Orange is for identifying the procedure.

**1 color for practitioner:**  
Vibrant Blue shows ownership to the practitioner.

## Did you know?


*The diagonal setup of the orange Rings shows the order of how the instruments should be used. This eliminates confusion between the doctor and assistant.*

Place your instrument to determine which ring size you need.



## E-Z ID Rings

- Prepackaged in an easy-to-use dispenser
- Available in Small, Large, and XL sizes to fit all instrument handles



### 25pk

Available in all colors.

- SM (70Z100) **\$7.50**
- LG (70Z200) **\$7.75**
- XL (70Z000) **\$8**

### Jewel System 200pk

Assorted colors (E, J, T, I, M, and U).

- SM (70Z115) **\$50**
- LG (70Z215) **\$51**
- XL (70Z015) **\$52**

### Vibrant System 200pk

Assorted colors (N, O, P, Q, R, and S).

- SM (70Z110) **\$50**
- LG (70Z210) **\$51**
- XL (70Z010) **\$52**

### Classic System 200pk

Assorted colors (A, B, D, G, and I).

- SM (70Z105) **\$50**
- LG (70Z205) **\$51**
- XL (70Z005) **\$52**



## E-Z ID Tape

- Tape dispenser has built-in tape cutter.
- Dispenser roller prevents tape from twisting.



### E-Z ID Tape Roll

Available in all colors.  
Includes : 1 - 10' roll.

**\$10.50** (70Z300)



### E-Z ID Tape System - Jewel

Assorted colors (E, J, T, I, M, and U).  
Includes : 8 - 3' rolls.

**\$43** (70Z315)

### E-Z ID Tape System - Vibrant

Assorted colors (N, O, P, Q, R, and S).  
Includes : 8 - 3' rolls.

**\$43** (70Z310)

### E-Z ID Tape System - Classic

Assorted colors (A, B, D, G, and I).  
Includes : 8 - 3' rolls.

**\$43** (70Z305)

Label Containers!  
Keep containers organized  
by labeling them per content.





# ISOLATION360™ POWER PACK

Zirc has the most extensive and innovative line of isolation products in the industry. There are practices that may be using one isolation device, not knowing there is something else that could be more beneficial to the procedure. With the Isolation360™ Power Pack, practices are able to explore and utilize different isolation products finding the best fit for their needs. We have an isolation device for every situation!





## Isolation360™ Power Pack

Includes:

- 1 Ideal Usage Chart
- 5 Lingua-Fix™
- 20 SE Cushion™
- 10 Pink Petal®
- 9 Airway Armor® (3 Small, 3 Medium, 3 Large)
- 10 Mr. Thirsty® One-Step (5 Peto/Small & 5 Med/Large)
- 10 Insti-Dam® (5 Original & 5 Relaxed Fit)

**\$100 (50Z950)**











































See pages 22-27 for individual product descriptions


**ISOLATION360™**  
POWER PACK

## WHAT PRODUCT SHOULD I USE?

\*Product image shows ideal option  
○ = Suggested Alternative

	RESTORATIVE						ENDODONTICS		ORAL SURGERY		HYGIENE		ORTHODONTICS					
	Crown & Bridge Inlays/Onlays	Crown & Bridge Inlays/Onlays	Composite / Amalgam	Digital Imaging	Veneer Prep	Veneer Insert	Anterior Tooth	Posterior Tooth	Extractions	Implant	Prophy	Scaling & Root Planing	Sealants	Ultrasonic Scaling	Bonding Single Bracket	Full Mouth Placement	Removal	Whitening (In-Office)
						 	 											
						 												
							<small>*Use in conjunction with Insti-Dam®</small>											
																		
						 												
	 <small>*for cementing</small>																	
																		

To see the full usage chart, visit...  
[www.zirc.com/360](http://www.zirc.com/360)



## A small solution to a BIG problem!

Dental items have been reported as the second most commonly ingested/aspirated foreign objects by adults.

Photo by:  
Jesus Escareno





Aspiration & Ingestion Prevention | Patient Protection

There are a multitude of small items used in dentistry: burs, files, implants and veneers to name a few. Combine these items with the moist environment of the oral cavity and you have a realistic probability of something being aspirated or ingested. Airway Armor® was created to eliminate the unnecessary risk for patients and protection for the practice.

**Ideal for Crown & Bridge, Endodontics, Extractions, Implants, Inlays/Onlays, Orthodontics, Veneers**

Aspiration & Ingestion prevention.  
Patient protection & comfort.  
Easy placement & removal.  
3 convenient sizes - Small, Medium & Large.  
Single Use.  
Latex Free.



### 15pk Trail Kit

Includes 5 of each size:

- (5) Small
- (5) Medium
- (5) Large

**\$45** (50Z960TRI)

### Airway Armor® 25pk

**\$65** Small (50Z960SM)

**\$65** Medium (50Z960MD)

**\$65** Large (50Z960LG)

To see an in-use video, visit...

[www.zirc.com/airwayarmor](http://www.zirc.com/airwayarmor)



## Let your Saliva Ejector Blossom!




Created by hygienist, Karrie Velky, the Pink Petal® is a great isolation product for all dental practices looking to fulfill the need of hands-free isolation. Easily slide the Pink Petal® on any saliva ejector creating a device that will provide continuous suction. This product is great relief for patients who want to hold the suction themselves or choose not to swallow during a procedure.

### Ideal for Scaling & Root Planing, Ultrasonic Scaling, Sealants

- Holds suction tip in place.
- Protects & retracts cheek.
- Patient can close around suction with ease.
- Easily slides onto saliva ejector.
- Soft, comfortable material.
- Single Use.
- Latex Free.



**Pink Petal® 50pk**  
\$25 (50Z975)



To see an in-use video, visit...  
[www.zirc.com/pinkpetal](http://www.zirc.com/pinkpetal)





## The Efficiency of Hands-Free Isolation

Hands-free retraction, isolation and evacuation provided by one simple device! Attach Mr.Thirsty® to your existing HVE, perform the procedure, and dispose when done! Rather than time spent on sterilizing or tracking number of uses, it can be spent on additional procedures, daily tasks or catching a break! According to a time and motion study, Mr.Thirsty® can reduce chair time up to 29%!

- Built-in bite block for patient comfort
- Connects directly to HVE Valve
- Keeps work environment dry
- Fully Assembled
- Latex-Free



Pedo/Small Adult

WITH ONE STEP, ACHIEVE  
**SIMPLICITY**

### AFFORDABLE

no upfront investment

### DISPOSABLE

no time wasted in maintenance  
or sterilization

### ADAPTABLE

connects to existing HVE, allowing it to be  
used in any operatory

### COMFORTABLE

provides tongue and cheek retraction and  
includes a built-in bite block, creating  
comfort for you and your patient

### TRIMMABLE

available in two sizes that can be trimmed  
to fit, no need to overstock with a variety  
of sizes



mr. thirsty® one-step

REDUCE CHAIR TIME  
BY AS MUCH AS  
**29%**



#### 10pk Trial Kit

Includes 5 of each size:

- (5) Pedo/Small Adult
- (5) Medium/Large

**\$35** (50Z987)



#### 25pk

- Pedo/Small Adult  
**\$70** (50Z986SM)
- Medium/Large  
**\$70** (50Z986LG)

ALSO  
AVAILABLE IN  
**100 PACK**

#### 100pk

- Pedo/Small Adult  
**\$225** (50Z986SM-B)
- Medium/Large  
**\$225** (50Z986LG-B)

Take comfort to  
the next level with  
lightweight tubing!



#### Comfort Kit

The Comfort Kit's lightweight tubing allows your heavy HVE Valve/Tubing to stay mounted to the equipment.

Includes:


- Lightweight lever HVE Valve
- 6' lightweight tubing
- mr. thirsty® one-step NOT included

**\$65** (50Z992)



Replacement HVE  
Valve can be found  
on page 37






To read a full product evaluation visit...  
<http://www.dentalproductshopper.com/relaxed-fit-insti-dam>

**Easy to use, non-threatening,  
 and comfortable for patient.**

## Insti-Dam®

### Features for BOTH Relaxed Fit and Original Insti-Dam®

- Built-in flexible medium gauge material that comes pre-assembled with frame and pre-punched hole
- Pre-punched hole makes positioning easy and allows eyes and nose to be open and comfortable
- Helps eliminate tearing (additional holes can be punched)
- Compact design fits outside patients' lips, easy-to-use, non-threatening and comfortable for patients
- Radiographs may be taken by folding it to the side
- Single use



### Latex-Free Insti-Dam® Relaxed Fit 20pk

Ideal for posterior isolation.  
 The extra "relaxed" material in the Relaxed Fit Insti-Dam® makes for quick posterior placement - no stretching, tearing, or pulling.  
 Dimensions: 4-7/8" x 4-7/8" x 6-1/8" diameter

**Patent #7044738**

**\$70 (50Z457)**

### Blue Latex-Free Insti-Dam® 20pk

Ideal for anterior isolation.  
 Made with polyisoprene.  
 Dimensions: 4-1/4" diameter

**\$60 (50Z459)**

### Natural Latex Insti-Dam® 20pk

Ideal for anterior isolation.  
 Dimensions: 4-1/4" diameter

**\$55 (50Z455)**



## Lingua-Fix® & SE Cushion



### Lingua-Fix®

Lingua-Fix® provides a dry work area.  
Holds tongue in a steady position to protect from injuries.  
Simple operation...connects to low volume suction.  
Comfortable, no sharp edges or corners.  
Large suction surface avoids blockage and is more comfortable.  
Reduces sedimentation of drilling by-products in pipes and tubing.  
Single use.  
Comes pre-assembled.

**\$13.50 (50Z953S) 10pk**

**\$49 (50Z953) 50pk**



### SE Cushion

Used in conjunction with the Saliva Ejector.  
Prevents tissue grabbing and tissue irritation.  
Increases patient comfort.  
Improves range of suction.  
Quiets the noise of the suction.  
Single use.

**\$22 (50Z970) 100pk**



**HELPFUL TIP:**  
Ideal for use  
with rubber  
dams

**For comfortable and efficient  
suction or evacuation.**

### The Relax Factor:

The extra "relaxed" material in the Relaxed Fit Insti-Dam® makes for quick posterior placement - no stretching, tearing or pulling.



### Insti-Dam® Dispenser

#### Made for Insti-Dam® Original ONLY

Can be mounted to wall or cabinet or can sit on a counter.  
Available in White (A) for latex or Neon Blue (N) for latex-free.  
Non-slip rubber bottom for stability.  
Holds 35 Insti-Dams®.  
Dimensions: 4-7/8" x 4-7/8" x 6-1/8"

**\$20 (50Z471) Color Code: A & N**

AM S



# High Definition Imagery




Experience high-definition quality. The brighter the workspace, the better the illumination and visibility. With ZIRC's Crystal HD® lens, you are visually working in a 40% brighter environment compared to the Rhodium lens that is currently the industry standard. If you haven't made the change to HD - what are you waiting for?

**40%  
BRIGHTER**  
Than Rhodium  
Mirrors



## CUTTING EDGE IMAGERY



**Rhodium = Reflects color 75% accurately**  
**Crystal HD® = Reflects color 99% accurately**

**Provides exceptionally bright,  
shadow-free images!**





Crystal HD®

Available in multiple styles!

Designed to reduce eye strain and fatigue.

Autoclavable

”

*As a microscope user, I can understand the importance of clear vision. ZIRC mirrors significantly improved the brightness of my image, giving my eyes highest quality vision. It's true - brightness makes all the difference.*

- Dr. Dimitris Tsanaksidis, Microscopic & Aesthetic Dentistry

## CRYSTAL HD® MIRROR + MIRROR MAGIC® The Perfect Pair

Pair our Crystal HD® lens with our mirror magic® System for the highest quality vision and best procedural outcome. See page 32.

CRYSTALHD®

+   
mirror magic®  
anti-fog system





# Crystal HD® Mouth Mirrors & Handles

Flush surface, no debris build-up can occur

Resin

Gap where debris build-up occurs

Stainless Steel

	RESIN	STAINLESS
Autoclavable	X	X
Antimicrobial protection	X	
Lens flush with mirror frame	X	
No galvanic shock on patient	X	
Rust-free	X	

**Why use  
Zirc's Resin  
Mirrors?**



## Assorted Soft Grip 12pk

### Vibrant Assorted

2 of each color: N, O, P, Q, R, and S.

**\$76 - #4** (50Z363-NEON)

**\$76 - #5** (50Z364-NEON)

### Jewel Assorted

2 of each color: E, J, T, I, M, and U.

**\$76 - #4** (50Z363-JWL)

**\$76 - #5** (50Z364-JWL)



## Assorted Thin Grip 12pk

### Vibrant Assorted

2 of each color: N, O, P, Q, R, and S.

**\$76 - #4** (50Z353-NEON)

**\$76 - #5** (50Z354-NEON)

### Jewel Assorted

2 of each color: E, J, T, I, M, and U.

**\$76 - #4** (50Z353-JWL)

**\$76 - #5** (50Z354-JWL)



## Black/Gray Soft Grip 12pk

**\$76 - #4** (50Z363)

**\$76 - #5** (50Z364)

## Black/Gray Soft Grip 4pk

**\$28.50 - #4** (50Z361)

**\$28.50 - #5** (50Z362)



## Handles



## Soft Grip Handle

**\$16 ea** (50Z340) Available in all colors




## Ergo Grip Handle

**\$16 ea** (50Z342) Available in all colors and K (black)



## What's your mirror size?



**#0** - 14mm

**#3** - 20mm

**#4** - 22mm

**#5** - 24mm

*diagram is to scale*



## Resin Screw-on Heads 12pk

Flush surface, no debris build-up can occur.

### 1-Color 12pk

Available in all colors and K (black).

**\$76 - #4** (50Z371)

**\$76 - #5** (50Z372)

### Vibrant Assorted

2 of each color: N, O, P, Q, R, and S.

**\$76 - #4** (50Z371-NEON)

**\$76 - #5** (50Z372-NEON)

### Jewel Assorted

2 of each color: E, J, T, I, M, and U.

**\$76 - #4** (50Z371-JWL)

**\$76 - #5** (50Z372-JWL)

### Classic Assorted

A, B, D, G, I, and K (black).

**\$76 - #4** (50Z371-PSTL)

**\$76 - #5** (50Z372-PSTL)



## Vibrant Assorted Double-Sided Mirror Head 4pk

1 of each color: N, O, R, and S.

**\$40 - #4** (50Z374-NEON)



## Stainless Steel Head 12pk

**\$76 - #3** (50R379)

**\$76 - #4** (50R380)

**\$76 - #5** (50R381)



## Gray Double-Sided Mirror Head 4pk

**\$40 - #4** (50Z374)



## Stainless Steel Head #0

**\$12.50 ea - #0** (50R378)



**ZIRC's Mirrors & Handles are autoclavable!**

## mirror magic® Anti-Fog System



without mirror magic®



with mirror magic®

**m**  
**w**  
mirror magic®  
anti-fog system

### Increase Your Chairside Efficiency

The mirror magic® Anti-Fog System prevents fogging on mouth mirrors and offers the practitioner a quick method for removing debris during procedures. The disposable Swipe Pads stick to the back of a gloved hand or a patient bib for easy access. The Anti-Fog System requires only a small amount of solution and is economic in size, making it easy to store.



Watch the magic in action  
at [www.zirc.com/mirrormagic](http://www.zirc.com/mirrormagic)



### mirror magic® Trial Kit

2-6ml bottles of Anti-Fog Solution  
24 Swipe Pads

\$25 (50R200)

Combine mirror magic® with Crystal HD®  
Mouth Mirror for the best view in dentistry.





# mirror magic® Anti-Fog System



*"The mirror magic® Anti-Fog System decreases fogging while also decreasing the build-up of debris on my mirror. This affords me outstanding visualization, which is crucial to performing quality dentistry. mirror magic® has truly made a difference in my practice and is one of those features I didn't realize I needed; and now I can't imagine practicing without it."*

— Andrew M. Goldsmith, DDS, DICOI, FIALD



## 10pk Anti-Fog Solution (6ml ea.)

1 bottle lasts approx. 50 patients.

**\$75** (50R205)



## Swipe Pads 96pk

Single use.

**\$25** (50R210)



## mirror magic® Organizer

Holds: 48 Swipe Pads and 1 - 6ml bottle of Anti-Fog Solution.

Store on operatory counter or in drawer.

Non-slip bottom for secure placement.

Note: box only, does not include mirror magic.

Dimensions: 3-3/4" x 3-5/8" x 1-7/8"

**\$16** (50R201)

S



**\$122 Value!**

**m**  
**w**  
mirror magic®  
anti-fog system



**CRYSTALHD®**

## Clarity Collection

Includes:

Anti-Fog Solution 10pk

Swipe Pads 96pk

1 Crystal HD® Mouth Mirror

1 mirror magic® Organizer

**\$100** (50R220)

## The Perfect Pair

The superior visual acuity and 40% brighter Crystal HD® lens combined with the quick debris removal and Anti-Fog Solution of mirror magic® have proven to be the perfect pair. Increase chairside efficiency and speed up your process with this dynamic duo.





## CAD CAM BLOCK LOCKER™

### Features

- Stores and manages CAD/CAM milling blocks and accessories
- Drawer accommodates accessories such as burs & lubricants
- Accommodates all blocks
- Wall mounted unit features a dark gray powder-coated finish for a modern look
- Provides a professional platform to present and discuss same-day restorations
- Online customizable tags to create organization by color and shade (<https://www.zirc.com/products/cad-cam-block-locker>)



## CAD/CAM Block Locker™ Large

Keep blocks in original packaging.  
Holds up to 440 blocks in a combination of packages and individual blocks.

Dimensions: 25" W x 26-5/8" T x 9-1/4" D

Weight: 49.4lbs

**\$975** (50Z416) Color Code: E, J, N, Q, S, T, U

**CONTROL BLOCK CHAOS**  
SHOWCASE YOUR MILLING BLOCKS IN STYLE



## CAD/CAM Block Locker™ Small

Designed to store individual/loose blocks only.  
Holds up to 288 individual blocks.

Dimensions: 19-1/2" W x 16" T x 9-1/4" D

Weight: 28.73lbs

**\$750** (50Z415) Color Code: E, J, N, Q, S, T, U

Pick Your Color!



E



J



N



Q



S



T



U



## Denture/Retainer Boxes



### Retainer Box

One piece construction.  
Solid locking mechanism.

#### Retainer Box 1" Deep 12 pk

Dimensions: 3-1/8" x 2-3/4" x 1"

**\$13.25 (25R500) 12pk**

Color Code: J, N, O, P, R, S, Assorted

#### Retainer Box 1-1/2" Deep 12 pk

Dimensions: 3-1/8" x 2-3/4" x 1-1/2"

**\$14.25 (25R550) 12pk**

Color Code: J, N, O, P, R, S, Assorted

S

High quality full color  
personalized imprinting  
now available!



### Denture Box

One piece construction.  
Solid locking mechanism.

#### Denture Box 12 pk

Dimensions: 4-3/8" x 3-3/8" x 2"

**\$14.75 (30R800) 12pk**

Color Code: E, I, J, M, T, U, Assorted

S

## NEW! Personalized Imprinting

Our NEW and IMPROVED imprinting service is quick and easy. Just select a full color icon from our Fun or Professional sets or upload your own logo. Add four lines of text for your business name, phone number, website address, or your own unique message. It's as easy as that!

### Upload Your Own Logo Option

Actual size of imprinting area



Brand  
your practice!

### Denture Box Imprinting 24 pk

Dimensions: 4-3/8" x 3-3/8" x 2"

**\$41 (30P800) 24pk**

Color Code: A, E, I, J, M, T, U, Jewel Assorted

### Retainer Box 1" Deep 24 pk

Dimensions: 3-1/8" x 2-3/4" x 1"

**\$38 (25P500) 24pk**

Color Code: A, J, N, O, P, R, S, Vibrant Assorted

### Retainer Box 1-1/2" Deep 24 pk

Dimensions: 3-1/8" x 2-3/4" x 1-1/2"

**\$40 (25P550) 24pk**


Color Code: A, J, N, O, P, R, S, Vibrant Assorted



### Fun Imprinting Options



### Professional Imprinting Options



Visit us online to order!

## Saliva Ejector Accessories & MEGA Products



HVE Valve - Lever



Saliva Ejector Valve - Lever

**Hands free isolation!**  
**See page 20-27.**

### Saliva Ejector Valve (with Shut-Off)



Quick disconnect with swivel.  
Rubber tip with protective screen.  
Provides variable control from full on to full off.

### Saliva Ejector Valve with Lever On/Off Control

Dimensions: 3" x 5/8" x 3/4"

**\$35 (23E365) Gray Only**

### HVE Valve (Lever style)

Dimensions: 3" x 5/8" x 11/16"

**\$45 (23E300) Silver Only**

### Saliva Ejector Valve Lever Only

Dimensions: 5/8" x 3/16" x 5/8"

**\$10 (23E364) Gray Only**

### Saliva Ejector Replacement Tip

Dimensions: 5/8" x 5/8" x 1"

**\$9 (23E363) Gray Only**



### Saliva Ejector Screens



Single use.

### Saliva Ejector Screens

Dimensions: 3/8" Dia. x 3/8"

**\$10.50 (50Z952S) 25pk**

**\$72 (50Z952P) 250pk**

## MEGA Products

### MEGA Gracey

15-1/2" in length.  
Great teaching tool.  
A gift for any dental personnel.

**\$20 (50Z312) Mustard Only**



### MEGA Mirror

(Patient hand mirror).  
2-1/2" lens and 15-1/2" in length.  
Keeps children entertained.  
Great gift for any dental personnel.

**\$20 (50Z309) Color Code: N, P, Q, R, S**



**Fun gift ideas for any dental personnel. Keeps children entertained.**



# Water Wise & Steri-Soaker



## Water Wise Starter Kit

Holds 1 liter.  
Includes bottle, adaptor, lid, neon yellow and neon blue bubbles (water level indicator), and tube.  
**\$33** (23Z100)



## Water Wise Replacement Bottle

Holds 1 liter.  
Includes bottle, blue lid, and neon yellow and neon blue bubbles (water level indicator).  
**\$23** (23Z105)



## Water Wise vs. Traditional Water Bottle

	Water Wise	Water Bottle
Easy to clean	X	
Antimicrobial	X	
Easy-to-fill large mouth	X	
Quick release	X	
Level indicator color coded for: Sterile and distilled water, Medicament	X	
Strong	X	X
Easy to grip	X	
Sealable lid	X	
No threads to strip	X	
Pressure relief valve	X	

We strive to make things easier around the dental office.

## The Steri-Soaker

An efficient, high-volume germicide tray designed for maximum use of chemicals to clean/disinfect instruments. It is designed for easy cleaning and is large enough to accommodate ZIRC Steri-Containers, Cassettes, Bur Blocks, and Bur Guards.

### Steri-Soaker

Calendar on lid conveniently reminds staff when to change the solution; use color code rings to mark dates. Transfer basket's ribbed and perforated design allows quick drain-off and easy pick-up of instruments. Holds 3.78L of liquid; fill lines shown in quarts and liters. Hinged fitted cover. Safe-to-use chemicals: cleaners (enzymes, detergents), disinfectants (phenols)/cold sterilants for glutaraldehyde and chlorines. Dimensions: 13-7/8" x 9-5/8" x 5-1/8"

**\$92** (50Z805) **Color Code: B, I**



Designed for easy cleaning!



# E-Z Access Shelf/Accessories

## E-Z Access Shelf

E-Z mounting to standard equipment post (1-3/4" - 2-1/4"); No tools needed.  
Will not damage post's finish.  
E-Z adjustment to multiple positions.  
Contoured edge to contain spills and materials.  
Stain resistant finish; surface disinfect only.  
Weight capacity of 25 lbs.  
Workable tray area: 13-3/4" L x 10-1/4" W.  
Standard E-Z Access Shelf as is accommodates 1-3/4" to 2-1/4" post.  
Rotates 360°.

**\$125 (20Z420) Color Code: A, G, I**



**Provides convenient work space for auxiliary items.**



## E-Z Access Shelf Adapters

### E-Z Access Shelf Adapter for Smaller Posts

Adapts to posts 1-3/8" to 2-3/4".  
Dimensions: 2-1/8" x 1-7/8" x 3/4"

**\$5 (20Z422) Color Code: A, G, I**

### 3-1/2" Bolt Adapter Kit for Larger Posts

Accommodates post sizes from 2-1/4" to 2-3/8".  
Dimensions: 4" x 2" x 2"

**\$10 (20Z437) Color Code: A, G, I**

### 4-1/2" Bolt Adapter Kit for Larger Posts

Accommodates post sizes from 2-3/8" to 2-3/4".  
Dimensions: 4-1/2" x 2" x 2"

**\$10 (20Z438) Color Code: A, G, I**



## Cup Holder-Snap On

Fits all standard cups.  
No tools needed.  
Dimensions: 3-1/8" x 2-1/2" x 1"

**\$11 (20Z436I) Gray Only**



## E-Z Access Shelf Surface Mount

Mount on horizontal flat surface only.  
Dimensions: 3-3/4" x 3-3/4" x 2-7/8"

**\$17 (20Z424) Color Code: A, G, I**



**Our products are able to fit every need!**



## 8" Extension Arm

Enables shelf to extend over patient.  
Shelf has strength to hold 25 lbs. when connected to arm.  
May be used with Wall Mount, Surface Mount, or directly to post.  
Dimensions: 10" x 4" x 4-1/2"

**\$31 (20Z428) Color Code: A, G, I**



## Wall or Cabinet Mount

Allows you to mount the E-Z Access Shelf onto a wall or the side of a cabinet. 90° angle of mount permits 180° movement of shelf on cabinet or next to wall.  
Dimensions: 6-1/2" x 5-3/4" x 4-7/8"

**\$31 (20Z426) Color Code: A, G, I**



# Cabinet Trays & Lab Pans

AM S 

## Cabinet Trays



### No 16A – Hand Instruments (Regular Well)

Dimensions: 7-7/8" x 3-3/4" x 15/16"

\$8.50 (20Z201) Color Code: A, B



### No 16 – Hand Instruments (Deep Well)

Dimensions: 7-3/4" x 3-3/4" x 1"


\$8.50 (20Z202) Color Code: A, B



### No 17 – Clamps, Retainers, etc.

Dimensions: 7-7/8" x 3-3/4" x 3/4"

\$8.50 (20Z203) Color Code: A, B



### No 18 – Finishing Strips

Dimensions: 7-7/8" x 3-3/4" x 3/4"

\$8.50 (20Z204) Color Code: A, B



### No 19 – Pliers, Scissors, etc.

Dimensions: 7-7/8" x 3-3/4" x 3/4"

\$8.50 (20Z205) Color Code: A, B



### No 20 – Discs, Wheels, Stones

Dimensions: 7-7/8" x 3-3/4" x 3/4"

\$8.50 (20Z206) Color Code: A, B



Wall-Hanging



AM S 

## Lab Pans

Organizes and transports impressions, models, dentures, and crown and bridge appliances. Stainless-steel spring clip holder. Color-coded.

### Lab Pan

Dimensions: 7-1/2" x 5-3/4" x 7"

\$6.75 (17Z101) Color Code: A, B, D, G, M

### Wall-Hanging Lab Pan

Same as standard version but with strip on top back for hanging on any vertical surface.

Saves on installing expensive shelving.

Dimensions: 7-3/4" x 5-3/4" x 7"

\$6.75 (17Z102) Color Code: A, B, D, G, M



## Wall Strip Lab Pan

Drilled and countersunk with screws and wall anchors. 18" Wall Strip.

Holds 3 Lab Pans.


Dimensions: 18" x 7/8" x 3/8"

\$9.50 (17Z104)

S

Save on installing expensive shelving.

## Dispensers



### Deluxe Cotton Roll Dispenser

Heavy-duty no-skid base.  
Dispenses up to three rolls at a time.  
Extra sturdy construction.  
Dimensions: 4-3/16" x 2-1/16" x 4"

**\$32 (28R800) Color Code: A, B**

S



### Cotton Roll Holder

Stainless pins embedded in weighted bottom.  
Neatly holds cotton rolls in an upright position.  
Hinged lid helps keep cotton rolls clean, prevents spilling.  
Dimensions: 3-5/8" x 3-1/4" x 2-3/8"

**\$30 (28R900) Color Code: A, B**

S



### Pellet Dispenser

Spring loaded – pushes pellets to the surface.

#### Dispenser for Small Pellets (sizes 3, 4, 5)

Dimensions: 1-7/8" x 1-7/8" x 2-1/16"

**\$31 (28R825) White Only**

#### Dispenser for Medium Pellets (sizes 0, 1, 2)

Dimensions: 1-7/8" x 1-7/8" x 2-1/16"

**\$31 (28R850) White Only**

S

### Sponge Dispensers

Spring action releases pads one at a time.  
Bottoms are weighted to prevent tipping.

#### 2" x 2" Sponge Dispenser

Dimensions: 2-1/2" x 2-1/2" x 2-7/8"

**\$31 (28R950) Color Code: A, B**

#### 4" x 4" Sponge Dispenser

Dimensions: 4-1/2" x 4-1/2" x 2-7/8"

**\$35 (28R975) Color Code: A, B**

S



### Matrix Dispenser

One roll of 3/8" x 12 yards Matrix clear tape 0.002 gauge.  
Allows Matrix to be contoured (u-shaped) as it dispenses.  
Stainless steel cutter allows for desired length.  
Disposable.

Dimensions: 2-1/8" x 2-1/4" x 5/8"

**\$26 (55Z200) Color Code: A, B**

S



# E-Z Storage/Handi-Hopper

## E-Z Storage Tub System

Hinges allow cover to stay open on its own!  
Each Tub comes with one vertical divider and one half-depth adjustment divider.  
Hinged cover comes in white or clear.

### E-Z Storage Organizer (with white cover)

Dimensions: 5-7/8" x 4-7/8" x 5-3/4"

Rack sold separately.

**\$35** (29R605)

### E-Z Storage Organizer (with clear cover)

Dimensions: 5-7/8" x 4-7/8" x 5-3/4"

Rack sold separately.

**\$35** (29R610)

### Replacement Hinges (pair)

Dimensions: 7/8" x 1/2" x 1/8"

**\$5** (29R627)

AM S

Hinges allow cover to stay open

Rack holding angled Tub



## E-Z Storage Single Rack

Holds storage organizer Tub.  
Tub can be placed flat or at an angle for easy access of supplies.  
Wall or cabinet mounted with optional mounting bracket.  
Sturdy metal construction; plastic-coated.  
White, non-slip feet.  
Dimensions: 7" x 5" x 4-3/4"

**\$21** (29R650)

S

## E-Z Storage Mounting Bracket

White metal bracket for the E-Z Storage Tub single rack.  
Pre-drilled holes for cabinet or wall mount.  
Dimensions: 5-1/8" x 1" x 1/2"

**\$5** (29R651)

S



AM S

Mount this mini waste receptacle within arm's-reach of chairside for convenient disposal of contaminated waste.

## Handi-Hopper (with Mounting Bracket)

Small plastic waste receptacle for disposal of contaminated waste at chairside, keeping tray clean.  
Slotted mounting bracket (included) for easy removal and replacement.  
To use for infectious waste, order red (Biohazard) liners.  
To use for other waste, order clear liners.

### Handi-Hopper

Dimensions: 7" x 4" 6-3/8"

**\$19** (20Z404) Color Code: G, I

### Handi-Hopper Bracket

Dimensions: 2-1/4" x 1/4" x 1"

**\$3.50** (20Z408) Color Code: G, I



## Handi-Hopper Liners (Single-Use)

Handi-Hopper Liners 100/Box (Red Biohazard)

**\$23** (20Z405) 100pk

Handi-Hopper Liners 100/Box (Clear)

**\$23** (20Z410) 100pk

AM S

## Waste-Drop Ring

Waste-Drop Ring for cabinet top.  
Fits 5" hole.  
Made of self-lubricating, easy-to-clean resin.  
Dimensions: 5-5/8" Dia. x 3"

**\$14** (22Z100) Color Code: A, G, I



# Cuspidor/Dispos-A-Screen Traps & Needle Capper



**Reduce possibilities of infection and contamination.**

## Cuspidor Screens

Single use.

### Cuspidor Screens Large Universal – 2-1/2" (gray)

Adec, Ampco, Belmont, Dab Atlante, Dansereau, DCI Den-Tal-Ez, Forest, Knight, Marus, MDT, Proma, Ritter, Schein.

**\$18 (51Z205S) Gray 25pk**

**\$73 (51Z205V) Gray 144pk**

### Cuspidor Screens Small – 1-5/8" (white)

Adec, Cascade and Performer cuspidor bowls.

**\$18 (51Z200S) White 25pk**

**\$73 (51Z200V) White 144pk**

## Dispos-A-Screen Traps

Avoid contact with potentially hazardous material.  
Improves performance of suction system.  
Color coded by size. Single use.

**\$18 (25pk)**

**\$72 (144pk)**

### Blue – 2-1/8"

Adec models with black plastic canisters, plus newer models by Ampco, Belmont, DCI, Forest, Marus MDT, Schein.

**25pk (51Z045S), 144pk (51Z045V)**

### Green – 1-7/8"

Ampco, Belmont, Biotec, Chayes, Dabi-Alante, Dansereau, DCI, Den-Tal-EZ, Dentech, Forest, Knight, Pelton & Crane, Proma and others.


**25pk (51Z123S), 144pk (51Z123V)**

### Gray – 2-3/4"

Midwest 210 Models.

**25pk (51Z070S), 144pk (51Z070V)**

*Left to right: Blue, Green, Gray, Yellow, Beige, Mauve and Coral*



### Yellow – 2-1/4"

Adec, model with metal canisters.

**25pk (51Z040S), 144pk (51Z040V)**

### Beige – 2-7/8"

Den-Tal-EZ models with metal canister.

**25pk (51Z060S), 144pk (51Z060V)**

### Mauve – 1-1/4"

Adec, Cascade, and Performer.

**25pk (51Z080S), 144pk (51Z080V)**

### Coral – 2-5/8"

Pelton and Crane Spirit Unit with Plastic Canister.

**25pk (51Z050S), 144pk (51Z050V)**



## Needle Capper

One-handed capping and recapping.

Syringe is for use; at convenient position in Needle Capper.

Base of device is secured to tabletop or tray with adhesive strip (included).

Needle is uncovered at same time syringe is withdrawn from Capper.

After use, syringe is replaced for recapping.

Needle Capper available in Standard and Long sizes.

Either size tube will fit standard base.

### Standard Needle Capper

(base, one tube, adhesive tape)

Dimensions: 2-3/16" x 1-1/4" x 1-1/4"

**\$26 (27R700) White Only**

### Long Needle Capper

(base, one tube, adhesive tape)

Dimensions: 2-1/2" x 1-1/4" x 1-3/4"

**\$26 (27R704) White Only**



# Mighty Mixer Bowls, Spatulas & Organizers



Single-use  
disposable mixing  
bowls.



## Mighty Mixer Bowls 36pk (White Only)

Disposable, eliminating cleanup time.  
Fits securely into mixing bowl without slippage.  
Saves wear on mixing bowl.  
Free-standing for use with or without mixing bowl.  
Capable for use with electric Alginate mixers.  
Single-Use.  
Dimensions: 4-7/8" Dia. x 3-1/8"

**\$21 (50Z505) 36pk**

## Mighty Mixer Bowls 36pk (Assorted)

Same features as the white Mighty Mixer Bowls.  
Dimensions: 4-7/8" Dia. x 3-1/8"  
9 of each color: J, N, P, and R  
Single-Use.

**\$21 (50Z505-ASTD) 36pk**

## Mighty Mixer Spatula 4pk (Assorted)

Same features as the white Mighty Mixer Spatula.  
Dimensions: 7-5/16" x 3/16" x 1-1/4"  
1 of each color: J, N, P, and R

**\$20 (50Z501-ASTD) 4pk**



## Mixing Sticks 100pk

Single use.  
Dimensions: 4-1/2" long  
**\$21 (50Z521) 100pk**



## Cement Spatula

White plastic, double-end spatula.  
For fillings, cement, and pastes.  
Self-lubricating plastic for easy cleaning.  
Dimensions: 6-3/8" x 1/4" x 3/8"

**\$6.50 ea. (50Z502)**



## Delrin Spatula 10pk

Flexible flat blade.  
Non-stick.  
Resistant to chemicals.  
Great for mixing cements, acrylics, etc.  
Dimensions: 4-5/8" x 3/16" Dia.

**\$10 (50Z523) 10pk**



## Mighty Mixer Spatula

Self-lubricating for easy cleaning.  
White plastic spatula.  
Broad, medium flex blade.  
Ideal for alginates and plaster.  
Dimensions: 7-5/16" x 3/16" x 1-1/4"

**\$6 ea. (50Z501) White Only**



## Implant Organizer

Quick and easy recognition. Individual compartments for implants.  
Shape of the organizer matches the arch of the mouth for implant placement.  
Separate compartments prevent implant appliances from mixing together.  
Lubricated plastic enables tooth number to be written in the designated area.  
Cover helps keep materials secure.  
Dimensions: 6-7/8" x 4-5/8" x 5/8"

**\$25 (20Z490B) Blue Only**



No more  
confusion!



## Restorative Organizer

Compartments shaped like the arch of the mouth for veneers, inlays, composite, and onlays.  
Disposable Mixing Wells included for bonding agent and cement.  
Accessory area for brushes and materials.  
Lubricated plastic enables tooth number to be written in the designated area.  
Light protective cover can be used as a light shield and to prevent materials from setting early.  
Includes 25pk of Mixing Wells.  
Dimensions: 6-7/8" x 4-5/8" x 5/8"

**\$30 (20Z495R) Neon Purple Only**



## Mixing Well (Single-Use) 144pk

Dimensions: 2-7/8" x 1-1/8" x 3/16"

**\$22 (20Z496V) 144pk**





MAKING DENTAL TEAMS  
HAPPIER & MORE *efficient!*



Because Time Is Everything

[www.zirc.com](http://www.zirc.com) | 800.328.2899 | [sales@zirc.com](mailto:sales@zirc.com)

\*All U.S. pricing subject to change